

Date .....

User name .....

Responsible .....

## Carbon fibre frame

### Wheelchair measures

Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48

The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.

Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP) ..... (Real distance between seat and footplate)

**IMPORTANT NOTICE!** From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)


Setting ☒ Standard ☐ Prudential ☐ Active ☐ Extreme

Backrest Height (TS) Std.Value 33,5 .....(min 24.5-max 47)

camber ☒ 0° ☐ 2° ☐ 4°

Front height STD 49 .....

Rear height STD 42 .....



☒ Version mudguard (Available only with rear wheels 24")

☐ Version straight sideguard

### Seat

☒ Seat cloth with pocket (under seat)

☐ Seat with tension straps 08883 £ 196

☐ Carbon rigid seat 08336 £ 415

☐ Pocket (under seat) [25x11x4,5] 08281 £ 65

### STYLE

☒ CARBON Essence


☐ ALA ROSSA 07947 £ 109

☐ ALA BIANCA 07948 £ 109

### Footplate

☒ Carbon fibre footplate - adjustable distance +/- 1,5cm. SPECIAL FLIP-UP!

Minimum footplate distance 35cm.



English  
Date ..... Order ..... User name ..... Address .....

**Backrest with rake adjustable from -4° to +16°**

- ☆ Backrest with double locking (upright and folded)  
☐ Backrest with upright locking system ONLY (no locking when folded)  
☐ Fixed back (non-folding back)

☆ Aluminium backrest

☐ Titanium backrest

☐ Carbon fibre backrest



06811 £ 109

(Lightest configuration)


04389 £ 600

- ☐ **A I R - Integrated posture backrest (not removable) - Lateral retaining 6,5cm.**  
 not available for seat width 48 - Std. Backrest upholstery not provided - Handle tube aluminium with short tubes

06815 £ 1015


**Lateral retaining wing B=low**

**Carbon shell H25**



Backrest Height (TS) min 29 ..... max 41


**Carbon shell H33**



Backrest Height (TS) ..... (37 + 46)

**Lateral retaining wing A=high**

**Carbon shell H30**



Backrest Height (TS) ..... (32 + 42,5)

**AIR Cover**

☆ Cover Standard - fabric insert: GREY, stitchings GREY

☐ **Cover SPORT**

09195 £ 87


- ☐ fabric insert: BLACK - Stitchings GREY  
☐ fabric insert: BLACK - Stitchings LIME  
☐ fabric insert: BLACK - Stitchings ORANGE  
☐ fabric insert: BLACK - Stitchings RED  
☐ fabric insert: BLACK - Stitchings BLUE  
☐ fabric insert: BLACK - Stitchings BLACK  
☐ fabric insert: BLACK - Stitchings FUCHSIA  
☐ fabric insert: BLACK - Stitchings ANTHRACITE

☐ **Cover ACTIVE**

09195 £ 87

- ☐ fabric insert: RED - Stitchings BLACK  
☐ fabric insert: LIME - Stitchings BLACK  
☐ fabric insert: BLUE - Stitchings BLACK

fabric insert:



Transpiring fabric stitchings

**Handle tube**

☆ Handle tube STD

- ☐ Handle tube with short tubes  
☐ Handle tube SPORT Without push handles

☐ Carbon fibre handle (with short tubes)

(Lightest configuration)

04390 £ 109

**Not available with: Posture backrest PHYSIO - Posture backrest A I R**

☐ Folding handles aluminium



02958 £ 218

**Casters**

☆ 4" Soft

☐ 3" Sport

(Lightest configuration)

☐ 5" Soft

☐ 3" Active Soft (wide) (Aluminium alloy rim)

06824 £ 109

☐ 4" Active Sport (narrow) (Aluminium alloy rim)

01398 £ 109

☐ 4" Active Soft (wide) (Aluminium alloy rim)

03943 £ 109

☐ 5" Active Soft (wide) (Aluminium alloy rim)

07518 £ 109

English

Date ..... Order ..... User name ..... Address .....

## Rear wheels

**Available with : VERSION MUDGUARD - VERSION STRAIGHT SIDEGUARD**

☆ 24" Ultra-light Progeo - Standard!

<input type="radio"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	(Lightest configuration)	04383	£ 742
<input type="radio"/> 24" Ultra-light "Spinergy LX"		01685	£ 742
<input type="radio"/> 24" Ultra-light "Spinergy Flex rim"		02270	£ 1200
<input type="radio"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes <input type="checkbox"/> Hub and blade spokes black	07846	£ 742
<input type="radio"/> 24" Off-road (aluminium handrim black)		01392	£ 164

**Available only with: VERSION STRAIGHT SIDEGUARD**

<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 26" Progeo - Standard!			
<input type="radio"/> 25" Ultra-light "Spinergy LX"		02101	£ 742
<input type="radio"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes	07847	£ 742
<input type="radio"/> 26" Ultra-light "Spinergy LX"		02102	£ 742

## Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	(Lightest configuration)		
<input type="checkbox"/> High pressure profile tyre	(24")			
<input type="radio"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669	£ 55
<input type="radio"/> Solid tyre Kik Mako - Colour Black	(24")		01041	£ 87
<input type="radio"/> Solid tyre Kik Mako - Colour Grey	(24")		01751	£ 87

## Handrim - Wheel distance

☆ Aluminium	☆ 3 cm	<input type="checkbox"/> 4 cm		
	(22"-24"-25"-26")	(Lightest configuration)		
<input type="radio"/> Aluminium black SOFT-TOUCH	(24"-25")		09170	£ 109
<input type="radio"/> Chromed Steel	(24"-25")		01684	£ 109
<input type="radio"/> Titanium	(24"-25")		00016	£ 218
<input type="radio"/> Ultra Grepp	(24"-25")		00017	£ 262
<input type="radio"/> Max Grepp	(24"-25")		00834	£ 305
<input type="radio"/> Ergo Para	(24")		00632	£ 305
<input type="radio"/> SURGE LT	(24"-25")	Oval countour with a rubber insert	06427	£ 316
<input type="radio"/> Quad (24")	(24")		00693	£ 175
<input type="radio"/> Silicon cover for standard aluminium handrim		Available for measures (22"-24"-25"-26")	04050	£ 142

## Wheel locks



☆ Without brake	(Lightest configuration)			
<input type="checkbox"/> Bent lever push to lock				
<input type="checkbox"/> Straight lever Push to lock				
<input type="checkbox"/> Straight lever Pull to lock				
<input type="radio"/> Sport brake OUT-FRONT			06428	£ 109
<input type="radio"/> Brake with extension lever push to lock			01949	£ 55
<input type="radio"/> Brake with extension lever pull to lock			03372	£ 55

## Armrests

<input type="radio"/> Swing away armrests	(Swing away, removable and height adjustable armrest)		00694	£ 218
	☆ Pad STANDARD			
	<input type="radio"/> Ergonomic pad (30x6 cm.)		04104	£ 55

Date .....	Order .....	English User name .....	Address .....
------------	-------------	----------------------------	---------------

**Optional 1**

- |   |  |   |                |                |
|---|--|---|----------------|----------------|
| <input type="radio"/> Transit wheels                |  |  | 03290          | £ 131          |
| <input type="radio"/> Anti-tipper with rapid hookup | (inside rotation to reduce the encumbrance)<br><input type="radio"/> Right<br><input type="radio"/> Left |  | 06821<br>06822 | £ 109<br>£ 109 |

**Optional 2**




- |  |  |  |       |       |
|--|--|--|-------|-------|
| <input type="radio"/> Pad for front frame protectors |  |  | 07058 | £ 39  |
| <input type="radio"/> Lumbar support                 | (to insert in the backrest upholstery) |  | 00013 | £ 33  |
| <input type="radio"/> Tetraclip                      | (Quick release axle for quadriplegic)  |  | 03947 | £ 131 |

**Additional Kit**

- |  |  |  |       |       |
|--|--|--|-------|-------|
| <input type="radio"/> 24" Off-road rear wheels | Off-Road tyre - aluminium handrim (black colour) |  | 03378 | £ 436 |
| <input type="radio"/> 24" FAT rear wheels      | FAT tyre 3" - aluminium handrim (black colour)   |  | 09357 | £ 764 |

**PHYSIO - Posture backrest**

Std. Backrest upholstery not provided

- |   |   |  |                |                |
|---|---|--|----------------|----------------|
| <input type="radio"/> Aluminum shell - Lateral retaining STANDARD 7,5cm               |    | Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41<br>Height <input type="checkbox"/> 52 | 05308<br>05309 | £ 884<br>£ 938 |
| Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low |   |  |                |                |
| <input type="radio"/> Aluminum shell - Lateral retaining DEEP 15cm                    |    | Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41<br>Height <input type="checkbox"/> 52 | 05310<br>05311 | £ 938<br>£ 993 |
| Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low |   |  |                |                |
| <input type="radio"/> Carbon shell - Lateral retaining STANDARD 7,5cm                 |  | Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41                                       | 05312          | £ 1015         |
| Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low |   |  |                |                |


**COVER**

☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




- ☐ SWING AWAY - Lateral support  
Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width,  
adjustable (0 ÷ -1,5cm)




- |  |                                      |
|--|--------------------------------------|
| Large (17,5x11,5)  |                                      |
| <input type="checkbox"/> Left 05646  | <input type="checkbox"/> Right 05647 |
| Small (14x9,5)   |                                      |
| <input type="checkbox"/> Left 05642  | <input type="checkbox"/> Right 05643 |
| <input type="checkbox"/> Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm) |                                      |

Internal thoracic support to the backrest width,  
adjustable (-4 ÷ -5,5cm)



- |                                     |                                      |
|-------------------------------------|--------------------------------------|
| Large (17,5x11,5)                   |                                      |
| <input type="checkbox"/> Left 05646 | <input type="checkbox"/> Right 05647 |
| Small (14x9,5)                      |                                      |
| <input type="checkbox"/> Left 05642 | <input type="checkbox"/> Right 05643 |

Internal thoracic support to the backrest width,  
adjustable (-8,5 ÷ -10cm)



- |                                     |                                      |                |
|-------------------------------------|--------------------------------------|----------------|
| Large (17,5x11,5)                   |                                      |                |
| <input type="checkbox"/> Left 05646 | <input type="checkbox"/> Right 05647 | £ 175<br>£ 175 |
| Small (14x9,5)                      |                                      |                |
| <input type="checkbox"/> Left 05642 | <input type="checkbox"/> Right 05643 | £ 175<br>£ 175 |
| (with additional kit)               |                                      |                |

- ☐ Headrest adjustable - Double ball mounting system  
Not available for Physio backrest's height 25 and with carbon shell

- ☐ Large (34x14)  
☐ Small (28x14)



05648 £ 393

English  
Date ..... Order ..... User name ..... Address .....

## REFERENCE TABLES

1	Rear wheels	
	Model	Features
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)

Italia Certified company normative UNI EN ISO 9001

Important notice! is page only if any optional have been chosen

Pag. 5 / 5

Wheelchair : CARBOMAX - Rev. 2018.3

Prices in EURO - VAT excluded  
to be added to the price list

☒ standard (supplied if no requests are indicated)  
☐ on request (without extra charge)  
☐ optional (cost to be added at the basic price)

Printed on 07-2018