

Date

User name

Company stamp

Responsible

Wheelchair measures

Seat width (LS) ☐33 ☐36 ☐39 ☐42 ☐45 ☐48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.

Seat depth (PS) ☐35 ☐37,5 ☐40 ☐42,5 ☐45 ☐47,5 ☐50
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP) (Real distance between seat and footplate)
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)


Setting ☆ Standard ☐ Prudential ☐ Active

Backrest Height (TS) Std.Value 33,5(min 24.5-max 47)

Camber ☆ 0° ☐ 2° ☐ 4°



Front frame

☆ Carbon fibre frame, Angle 90°



Front height - Tick the box with the required height

The possible difference in height between HIGH frame and LOW frame is 3cm (with same fork, caster and assembly hole).

	 SPORT fork (H95)	 Small size fork (H115)
LOW FRAME 90° Front height (min 44,5 - max 50) Footplate distance (min 34,5) Footplate distance (min 18) with U-TUBE Aluminium footplate	<p>Casters</p> <p>3" <input type="checkbox"/> ~44,5 <input type="checkbox"/> ~45,5</p> <p>4" <input type="checkbox"/> ~45,5 <input type="checkbox"/> ~46,5</p> <p>5" <input type="checkbox"/> ~48</p>	<p>Casters</p> <p>3" <input type="checkbox"/> ~46,5 <input type="checkbox"/> ~47,5</p> <p>4" <input type="checkbox"/> ~46,5 <input type="checkbox"/> ~47,5 <input type="checkbox"/> ~48,5</p> <p>5" <input type="checkbox"/> ~49 <input type="checkbox"/> ~50</p>
HIGH FRAME 90° Front height (min 47,5 - max 53) Footplate distance (min 37,5) Footplate distance (min 21) with U-TUBE Aluminium footplate	<p>Casters</p> <p>3" <input type="checkbox"/> ~47,5 <input type="checkbox"/> ~48,5</p> <p>4" <input type="checkbox"/> ~48,5 <input type="checkbox"/> ~49,5</p> <p>5" <input type="checkbox"/> ~51</p>	<p>Casters</p> <p>3" <input type="checkbox"/> ~49,5 <input type="checkbox"/> ~50,5</p> <p>4" <input type="checkbox"/> ~49,5 <input type="checkbox"/> ~50,5 <input type="checkbox"/> ~51,5</p> <p>5" <input type="checkbox"/> ~52 <input type="checkbox"/> ~53</p>

The real front-seat height could be slightly different from the chosen one considering it is affected by seat depth and rear seat height.

Rear height

The difference between Font Height and Rear Height must be in the 3 - 12 cm range.

Rear wheels 22" (min 36 - max 44)
Rear wheels 24" (min 38 - max 46)
Rear wheels 25" (min 39 - max 47)
Rear wheels 26" (min 41 - max 48)

Seat

☆ Seat cloth with pocket (under seat)

☐ Seat with tension straps 08883 £ 196

☐ Carbon rigid seat 08336 £ 415

☐ Pocket (under seat) [25x11x4,5] 08281 £ 65

Date Order English User name Address

CUSTOMIZE COLORS

① Fork colour

☆ Opaque black

☐ Bright silver

☐ Pearl white

☐ Bright red

08342

☐ Bright blue

☐ Bright fuchsia

☐ Bright lime

£ 55

② Colour rear wheels

To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)

☆ Black

☐ Bright silver

☐ Pearl white

☐ Bright red


04692

☐ Bright blue

☐ Bright fuchsia

☐ Bright lime

£ 55



Rear wheels axle

☆ Aluminum rear wheels axle

☐ Carbon fiber rear wheels axle

(Lightest configuration)

07647

£ 131

Backrest with rake adjustable from -4° to +16°

☆ Adjustable velcro straps and padded transpiring upholstery

☐ Light adjustable velcro straps and padded transpiring upholstery SPORT

08337

£ 65

☆ Backrest with double locking (upright and folded)

☐ Backrest with upright locking system ONLY (no locking when folded)

☐ Fixed back(non-folding back)

☆ Aluminium backrest

☐ Titanium backrest

☐ Carbon fibre backrest

(Lightest configuration)



06811

£ 109

04389

£ 600

☐ A I R - Integrated posture backrest (not removable) - Lateral retaining 6,5cm.


not available for seat width 48 - Std. Backrest upholstery not provided - Handle tube aluminium with short tubes

06815

£ 1015

Lateral retaining wing B=low

Carbon shell H25



Backrest Height (TS) min 29 max 41

Carbon shell H33



Backrest Height (TS) (37 ÷ 46)

Lateral retaining wing A=high

Carbon shell H30



Backrest Height (TS) (32 ÷ 42,5)

AIR Cover

☆ Cover Standard - fabric insert: GREY, stitchings GREY

☐ Cover SPORT

09195

£ 87

☐ fabric insert: BLACK - Stitchings GREY

☐ fabric insert: BLACK - Stitchings LIME

☐ fabric insert: BLACK - Stitchings ORANGE

☐ fabric insert: BLACK - Stitchings RED

☐ fabric insert: BLACK - Stitchings BLUE

☐ fabric insert: BLACK - Stitchings BLACK

☐ fabric insert: BLACK - Stitchings FUCHSIA

☐ fabric insert: BLACK - Stitchings ANTHRACITE

☐ Cover ACTIVE

09195


£ 87

☐ fabric insert: RED - Stitchings BLACK

☐ fabric insert: LIME - Stitchings BLACK

☐ fabric insert: BLUE - Stitchings BLACK

fabric insert:



Transpiring fabric stitchings

Date Order English User name Address

Handle tube

☆ Handle tube STD

- ☐ Handle tube SPORT Without push handles
☐ Handle tube with short tubes

☐ Carbon fibre handle (with short tubes) **Weight reduction approx. 110gr.**
Not available with: Posture backrest PHYSIO - Posture backrest A I R

04390 £ 109

☐ Folding handles aluminium




02958 £ 218

Footplate

☆ Aluminium footplate. Rake not adjustable

☆ Aluminium footplate ☐ with plastic plate


☐ One piece aluminium footplate - Flip-up! FLIP-BACK! - rake and depth adjustable



09420 £ 164

☐ Aluminium footplate, rake and depth adjustable


☆ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External



04456 £ 218

☐ U-TUBE - Aluminium footplate, rake and depth adjustable


Low Frame minimum footplate distance 18 - High Frame minimum footplate distance 21



09853 £ 196

☐ One piece Carbon fibre footplate. Rake adjustable

☆ Internal ☐ External



04110 £ 327

Mudguard

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard

☐ Removable version



08921 £ 109

Casters

☆ 4" Soft

- ☐ 3" Sport
☐ 5" Soft

☐ 3" Active Soft (wide) (Aluminium alloy rim) 06824 £ 109
☐ 4" Active Sport (narrow) (Aluminium alloy rim) 01398 £ 109
☐ 4" Active Soft (wide) (Aluminium alloy rim) 03943 £ 109
☐ 5" Active Soft (wide) (Aluminium alloy rim) 07518 £ 109

Rear wheels

☆ 24" Ultra-light Progeo - Standard!

- ☐ 22" Ultra-light Progeo - Standard!
☐ 25" Ultra-light Progeo - Standard! (colour: black)
☐ 26" Progeo - Standard! (colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT **Weight reduction approx. 450gr.**

04383 £ 742

☐ 24" Ultra-light "Spinergy LX"

01685 £ 742

☐ 24" Ultra-light "Spinergy LXL"

☐ Black hub and white blade spokes

☐ Hub and blade spokes black

07846 £ 742

☐ 24" Ultra-light "Spinergy Flex rim"

02270 £ 1200

☐ 25" Ultra-light "Spinergy LX"

02101 £ 742

☐ 25" Ultra-light "Spinergy LXL"

☐ Black hub and white blade spokes

07847 £ 742

☐ 26" Ultra-light "Spinergy LX"

02102 £ 742

☐ 24" Off-road (aluminium handrim black)

01392 £ 164

Rear wheel tyre

☆ High pressure tyre (22"-24"-25"-26")

☐ High pressure profile tyre (24")

☐ Marathon Plus (24"-25")

with anti-puncture protection

01669 £ 55

☐ Solid tyre Kik Mako - Colour Black (24")

01041 £ 87

☐ Solid tyre Kik Mako - Colour Grey (24")

01751 £ 87

English

Date Order User name Address

Handrim - Wheel distance

★ 3 cm □ 4 cm

☆ Aluminium	(22"-24"-25"-26")		
○ Aluminium black SOFT-TOUCH	(24"-25")	09170	£ 109
○ Chromed Steel	(24"-25")	01684	£ 109
○ Titanium	(24"-25")	00016	£ 218
○ Ultra Grepp	(24"-25")	00017	£ 262
○ Max Grepp	(24"-25")	00834	£ 305
○ Ergo Para	(24")	00632	£ 305
○ SURGE LT	(24"-25")	06427	£ 316
○ Quad (24")	(24")	00693	£ 175
○ Silicon cover for standard aluminium handrim		04050	£ 142

Oval countour with a rubber insert

Available for measures (22"-24"-25"-26")

Wheel locks

☆ Bent lever push to lock

□ Straight lever Push to lock

□ Straight lever Pull to lock

○ Sport brake OUT-FRONT

○ Brake with extension lever push to lock

○ Brake with extension lever pull to lock



06428 £ 109

01949 £ 55

03372 £ 55

Armrests

○ Swing away armrests

(Swing away, removable and height adjustable armrest)

○ Removable armrest with carbon-fiber straight sideguard



00694 £ 218

09079 £ 196

Kind of pad to be mounted on chosen armrest

☆ Pad STANDARD

○ Ergonomic pad (30x6 cm.)

04104 £ 55

Optional 1

Only one optional devices for side can be chosen

○ Transit wheels

Exclude the installation of Anti-tippers

○ Anti-tipper with rapid hookup

Exclude the installation of Transit wheels
(inside rotation to reduce the encumbrance)

○ Right

○ Left



03290 £ 131

06821 £ 109

06822 £ 109

Optional 2

○ Push handles for companion

(Removable and height adjustable)

○ Pad for front frame protectors

03645 £ 153

07058 £ 39

○ Tetraclip

(Quick release axle for quadriplegic)

03947 £ 131

○ Safety pelvic belt

03972 £ 55

○ Lumbar support

(to insert in the backrest upholstery)

00013 £ 33

○ Swing away lateral support

○ Right

○ Left

Not available with:

Posture backrest PHYSIO - Posture backrest A I R



02097 £ 218

02098 £ 218

○ Spoke guards Black

(22"-24")

01043 £ 131

○ Spoke guards Transparent

(22"-24")

00835 £ 131

○ Spoke guards Progeo design

(22"-24")

00019 £ 153

○ Spoke guards Dragon

(22"-24")

02106 £ 153

Additional Kit

○ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378 £ 436

○ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357 £ 764

Date Order English User name Address

PHYSIO - Posture backrest

☐ Aluminum shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41
Height ☐ 52

05308 £ 884
05309 £ 938



Lateral retaining wing ☐ A=high ☐ B=low

☐ Aluminum shell - Lateral retaining DEEP 15cm

Height ☐ 25 ☐ 33 ☐ 41

05310 £ 938

H25 not available for seat width 39-42-45-48



Height ☐ 52

05311 £ 993

Lateral retaining wing ☐ A=high ☐ B=low

☐ Carbon shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41

05312 £ 1015



Lateral retaining wing ☐ A=high ☐ B=low

COVER


☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems



☐ SWING AWAY - Lateral support

Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)



Large (17,5x11,5)


☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)




Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)



Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

(with additional kit)

£ 175
£ 175

£ 175
£ 175

☐ Headrest adjustable - Double ball mounting system

Not available for Physio backrest's height 25 and with carbon shell



☐ Large (34x14)
☐ Small (28x14)

05648 £ 393

English
Date Order User name Address

Rehateam s.r.l. Vicolo Negrelli, 5 - 31040 Castagnole di Paese (TV) - Tel. +39.0422.484657 / +39.0422.480387 - Fax +39.0422.484661 / +39.0422.484684
P.IVA IT 03182520266 - C.C.I.A.A. TV n. 225967 - Iscr. Trib. TV n. 41690 Capitale Sociale I.V. Euro 100.000,00

Email: info@rehateamprogeo.com - web: http://www.rehateamprogeo.com

1	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)

Note

Italia Certified company normative UNI EN ISO 9001

Pag. 6 / 6

Wheelchair : DUKE - Rev. 2019.0

Prices in EURO - VAT excluded
to be added to the price list

Important notice! s page only if any optional have been chosen

☒ standard (supplied if no requests are indicated)
☐ on request (without extra charge)
☐ optional (cost to be added at the basic price)