

Date

User name

Company stamp

Responsible

Wheelchair measures

Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.

Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP) (Real distance between seat and footplate)
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)

Setting ☒ Standard ☐ Prudential ☐ Active

Backrest Height (TS) Std.Value 33,5 (min 24.5-max 47)


camber ☒ 0° ☐ 2°

Front height STD 49


Rear height STD 42

Front frame

☒ Aluminium standard frame




☒



☐

(max 5" caster)



☐

(max 4" caster)

Colour aluminium frame

☒ Bright silver

- ☐ Polished black
☐ Opaque black
☐ Bright blue
☐ Bright anthracite
☐ Bright orange
☐ Bright fuchsia
☐ Bright red
☐ Bright lime
☐ Bright blue china
☐ Pearl white


☐ CUSTOM - Aluminium frame - welded made-to-measure

07652


£ 327

Available only with Footplate COMPACT, automatic Closure adjustable distance +/- 1,5cm (Lightest configuration)

Pay attention to all measures. EGO CUSTOM is a welded frame wheelchair, thus seat height and setting can not be modified after its manufacturing.




☐



☐

(max 5" caster)



☐

(max 4" caster)

Frame colour (without extra charge)

- ☐ Bright silver
☐ Opaque black

- ☐ Pearl white
☐ Polished black

- ☐ Opaque orange

Frame colour on request

- ☐ Metallic matt blue

- ☐ Metallic matt grey

- ☐ Metallic matt anthracite


09194

£ 327


☐ Carbon fibre frame

03347


£ 1309



☐ (max 5" caster)



☐ (max 5" caster)



☐ (max 4" caster)

☒ Polished transparent painting

☐ Opaque transparent painting

08922

£ 109

Backrest with rake adjustable from -4° to +12°

☒ With push handles short

☐ Sport Without push handles

☐ Backrest with folding handles

(Lightest configuration)



02958

£ 218

Date Order English User name Address

Casters

- ☆ 4" Soft
- ☐ 3" Sport (Lightest configuration)
- ☐ 5" Soft
- ☐ 6" Soft (Available only with 100° angle)
- ☐ 6" Pneumatic (Available only with 100° angle)
- ☐ 3" Active Soft (wide) (Aluminium alloy rim) 06824 £ 109
- ☐ 4" Active Sport (narrow) (Aluminium alloy rim) 01398 £ 109
- ☐ 4" Active Soft (wide) (Aluminium alloy rim) 03943 £ 109
- ☐ 5" Active Soft (wide) (Aluminium alloy rim) 07518 £ 109

Support front fork - Not available with ALUMINIUM FRAME CUSTOM

- ☆ Standard support front fork Exclude the installation of Front fork FROG LEGS



- ☐ Front fork support V Design



Special Options

- ☐ Rear wheels plate extension Not available with ALUMINIUM FRAME CUSTOM 01950 £ 131
- ☐ lockingband/strip to close/fold the wheelchair 06833 £ 22
- ☐ Drawstring to unlock the back with a singleaction 06834 £ 109

Colour set: Rear wheels and forks - To be matched only with:Rear wheels Ultra-light Progeo - Standard! (22 - 24)

- ☆ Black 22" 24" 25" 26"
- ☐ Bright silver 22" 24" 00833 £ 109
- ☐ Bright anthracite 24" Not available with ALUMINIUM FRAME CUSTOM 00833 £ 109
- ☐ Bright orange 22" 24" Not available with ALUMINIUM FRAME CUSTOM 00833 £ 109
- ☐ Bright fuchsia 22" 24" 00833 £ 109
- ☐ Bright blue 22" 24" 00833 £ 109
- ☐ Bright red 22" 24" 00833 £ 109
- ☐ Bright lime 22" 24" 00833 £ 109
- ☐ Pearl white 24" 00833 £ 109

Rear wheels

- ☆ 24" Ultra-light Progeo - Standard!
- ☐ 22" Ultra-light Progeo - Standard!
- ☐ 25" Ultra-light Progeo - Standard! (colour: black)
- ☐ 26" Progeo - Standard! (colour: black)
- ☐ 24" ULTRA-LIGHT PROGEO KEY SPORT (Lightest configuration) 04383 £ 742
- ☐ 24" Ultra-light "Spinergy LX" 01685 £ 742
- ☐ 24" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes ☐ Hub and blade spokes black 07846 £ 742
- ☐ 24" Ultra-light "Spinergy Flex rim" 02270 £ 1200
- ☐ 25" Ultra-light "Spinergy LX" 02101 £ 742
- ☐ 25" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes 07847 £ 742
- ☐ 26" Ultra-light "Spinergy LX" 02102 £ 742
- ☐ 24" Off-road (aluminium handrim black) 01392 £ 164

Rear wheel tyre

- ☆ High pressure tyre (22"-24"-25"-26") (Lightest configuration)
- ☐ High pressure profile tyre (24")
- ☐ Marathon Plus (24"-25") with anti-puncture protection 01669 £ 55
- ☐ Solid tyre Kik Mako - Colour Black (24") 01041 £ 87
- ☐ Solid tyre Kik Mako - Colour Grey (24") 01751 £ 87

Handrim - Wheel distance

- ☆ 3 cm ☐ 4 cm
- ☆ Aluminium (22"-24"-25"-26") (Lightest configuration)
- ☐ Aluminium black SOFT-TOUCH (24"-25") 09170 £ 109
- ☐ Chromed Steel (24"-25") 01684 £ 109
- ☐ Titanium (24"-25") 00016 £ 218
- ☐ Ultra Grepp (24"-25") 00017 £ 262
- ☐ Max Grepp (24"-25") 00834 £ 305
- ☐ Ergo Para (24") 00632 £ 305
- ☐ SURGE LT (24"-25") Oval countour with a rubber insert 06427 £ 316
- ☐ Quad (24") (24") 00693 £ 175
- ☐ Silicon cover for standard aluminium handrim Available for measures (22"-24"-25"-26") 04050 £ 142

Date Order English User name Address


Wheel locks

☆ Bent lever push to lock

- ☐ Straight lever Push to lock
☐ Straight lever Pull to lock

☐ Sport brake OUT-FRONT

(Lightest configuration)



06428

£ 109

Automatic Closure footplate

☆ Footplate COMPACT, automatic Closure

(Lightest configuration)

Minimum footplate distance 35cm.

RAKE AND DEPTH ADJUSTABLE

☐ Internal

☐ External



Footplate - See reference tables N°3 - Not available with ALUMINIUM FRAME CUSTOM


Warning!! All flip-up footrests, need a stretch calf band and standard fork-support

Weight increase approx. 500gr.

☐ Plastic footrest, rake adjustable

Minimum footplate distance 35cm.

☆2/3 Internal ☐ 2/3 External




☐ One piece footplate - Flip-up!, rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal

☆2/3 Internal

☐ External




☐ One piece aluminium footplate - Flip-up!, rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal

☆2/3 Internal

☐ 2/3 External



☐ One piece Carbon fibre footplate, rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.

☐ Internal

☐ External



03291

£ 240

Mudguard

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard

☐ Extra charge for : REMOVABLE VERSION



08921

£ 109

Removable armrest with ABS straight sideguard

☐ Removable armrest with carbon-fiber straight sideguard

☆ Pad STANDARD

☐ Ergonomic pad (30x6 cm.)



09079

£ 196

04104

£ 55

Optional 1

Important notice! With the EGO CUSTOM you may choose only one of these accessories with fee

☐ Swing away armrests

(Swing away, removable and height adjustable armrest)

☆ Pad STANDARD

☐ Ergonomic pad (30x6 cm.)



00694

£ 218

04104

£ 55

☐ Transit wheels



00836


£ 131

☐ Anti-tipper with rapid hookup

(inside rotation to reduce the encumbrance)

☐ Right

☐ Left



00023

£ 109

00024

£ 109

☐ Tipping aid

☐ Right

☐ Left



00025

£ 55

01046

£ 55

☐ Crutches support

☐ Right

☐ Left



00838

£ 65

00837

£ 65

Date Order English User name Address

Optional 2

<input type="checkbox"/> Tetraclip	(Quick release axle for quadriplegic)	03947	£ 131
<input type="checkbox"/> Safety pelvic belt		03972	£ 55
<input type="checkbox"/> Quick release front forks	Available with: Standard support front fork	00030	£ 131
<input type="checkbox"/> Pad for front frame protectors		07058	£ 39
<input type="checkbox"/> Push handles for companion	(Removable and height adjustable)	01048	£ 153
<input type="checkbox"/> Harness		04103	£ 109
<input type="checkbox"/> Lumbar support	(to insert in the backrest upholstery)	00013	£ 33
<input type="checkbox"/> Swing away lateral support	<input type="checkbox"/> Right	02097	£ 218
Not available with Posture backrest PHYSIO	<input type="checkbox"/> Left	02098	£ 218
<input type="checkbox"/> Spoke guards Black	(22"-24")	01043	£ 131
<input type="checkbox"/> Spoke guards Transparent	(22"-24")	00835	£ 131
<input type="checkbox"/> Spoke guards Progeo design	(22"-24")	00019	£ 153
<input type="checkbox"/> Spoke guards Dragon	(22"-24")	02106	£ 153






Additional Kit

<input type="checkbox"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	£ 436
<input type="checkbox"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	£ 764

PHYSIO - Posture backrest

Std. Backrest upholstery not provided

<input type="checkbox"/> Aluminum shell Lateral retaining STANDARD 7,5cm		Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41 Height <input type="checkbox"/> 52	05308 05309	£ 884 £ 938
		Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low		
<input type="checkbox"/> Aluminum shell Lateral retaining DEEP 15cm		Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41 H 25 not available for seat width 39-42-45-48 Height <input type="checkbox"/> 52	05310 05311	£ 938 £ 993
		Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low		
<input type="checkbox"/> Carbon shell Lateral retaining STANDARD 7,5cm		Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41	05312	£ 1015
		Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low		

COVER


☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems



SWING AWAY - Lateral support


Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)




<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647
		<input type="checkbox"/> Right	05647
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643
		<input type="checkbox"/> Right	05643
<input type="checkbox"/> Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)			

Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)



<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647
		<input type="checkbox"/> Right	05647
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643
		<input type="checkbox"/> Right	05643

Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)



<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647	£ 175
		<input type="checkbox"/> Right	05647	£ 175
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643	£ 175
		<input type="checkbox"/> Right	05643	£ 175
<input type="checkbox"/> Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)				

Date	Order	English User name	Address
------------	-------------	----------------------------	---------------

☐ Headrest adjustable - Double ball mounting system



Not available for Physio backrest's height 25 and with carbon shell

- ☐ Large (34x14)
☐ Small (28x14)



05648 £ 393

Note

Italia Certified company normative UNI EN ISO 9001			Important notices page only if any optional have been chosen	
Pag. 5 / 5	Wheelchair : EGO - Rev. 2019.0	Prices in EURO - VAT excluded to be added to the price list	<input checked="" type="checkbox"/> standard	(supplied if no requests are indicated)
			<input type="checkbox"/> on request	(without extra charge)
			<input type="radio"/> optional	(cost to be added at the basic price)

Enclosure to wheelchair order form Ego

REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SMALL 	max 46	max 47	max 48	
MEDIUM 	max 50	max 51	max 52	max 54
The size of the fork will be defined only by Rehateam Progeo in accordance with wheelchair final configuration.				

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 35	min 38	min 39	min 40	
max 43	max 45	max 46	max 48	

3	Footplate option	Footplate distance
	Footplate COMPACT, automatic Closure	min 36
	Plastic footrest	min 36
	One piece footplate - Flip-up!	min 36
	One piece aluminium footplate - Flip-up!	min 36
	One piece Carbon fibre footplate	min 36
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)		
NOTE! Lower footplate distances are possible only after our valuation		
Calculated with the best configuration and with the footplates positioned on inside.		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)