

Date .....

User name .....

Company stamp

Responsible .....

## Wheelchair measures

Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48

Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5

Footplate distance (DP) ..... **(Real distance between seat and footplate)**  
**IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)**

Setting ☆ Standard ☐ Prudential ☐ Active ☐ Extreme

Backrest Height (TS) Std.Value 36 .....(min 30-max 54)


Crossbrace (black) ☆ Single ☐ Double(Not seat width 33-36)

Front height ..... (STD 49) - Rear height ..... (STD 44)

## Front frame


- ☐ Short frame - Ideal with Seat depth (PS) 35 - 37,5  
☐ Medium frame - Ideal with Seat depth (PS) 40 - 42,5  
☐ Long frame - Ideal with Seat depth (PS) 45 - 47,5

## Frame Angle



## Rear frame

☆ Aluminium standard frame



camber

☆ 0° ☐ 1,5° ☐ 3°

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one.



○ Aluminium frame DYNAMIC

With rake adjustable backrest (±9°)

Black rear wheel plate with multifunction jointed system

camber

☆ 0° ☐ 1,5° ☐ 2° ☐ 4°

The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.


The real seat depth, from the front part of the seat cloth to the backrest tube will be 2 cm more than the declared one.

## Colour aluminium frame

- ☐ Polished black  
☐ Bright silver  
☐ Bright blue  
☐ Bright anthracite  
☐ Bright orange  
☐ Bright fuchsia  
☐ Bright red  
☐ Bright lime  
☐ Bright blue china  
☐ Pearl white

03942 £ 164

Aluminium frame  
Titanium frame  
TMAX



○ Titanium frame T MAX

**Brushed Titanium frame - not varnished**

**Weight reduction approx. 250gr**

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one.

03940 £ 327

## Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard!

☆ Polished black	22" 24" 25" 26"		
○ Bright silver	22" 24"	00833	£ 109
○ Bright anthracite	24"	00833	£ 109
○ Bright orange	22" 24"	00833	£ 109
○ Bright fuchsia	22" 24"	00833	£ 109
○ Bright blue	22" 24"	00833	£ 109
○ Bright red	22" 24"	00833	£ 109
○ Bright lime	22" 24"	00833	£ 109
○ Pearl white	24"	00833	£ 109

## Casters

☆ 5" Soft			
<input type="checkbox"/> 4" Soft			
<input type="checkbox"/> 3" Sport			
<input type="checkbox"/> 6" Soft	Not available with frame angle 95°		
<input type="checkbox"/> 6" Pneumatic	Not available with frame angle 95°		
<input type="checkbox"/> 7" Soft	Not available with frame angle 95°		
○ 3" Active Soft (wide)	(Aluminium alloy rim)	06824	£ 109
○ 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	£ 109
○ 4" Active Soft (wide)	(Aluminium alloy rim)	03943	£ 109
○ 5" Active Soft (wide)	(Aluminium alloy rim)	07518	£ 109

English

Date ..... Order ..... User name ..... Address .....

## Rear wheels

- ☆ 24" Ultra-light Progeo - Standard!
- ☐ 22" Ultra-light Progeo - Standard!
- ☐ 25" Ultra-light Progeo - Standard! (colour: black)
- ☐ 26" Progeo - Standard! (colour: black)
- 24" ULTRA-LIGHT PROGEO KEY SPORT **Weight reduction approx. 450gr.** 04383 **£ 742**
- 24" Ultra-light "Spinergy LX" 01685 **£ 742**
- 24" Ultra-light "Spinergy Flex rim" 02270 **£ 1200**
- 24" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes ☐ Hub and blade spokes black 07846 **£ 742**
- 25" Ultra-light "Spinergy LX" 02101 **£ 742**
- 25" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes 07847 **£ 742**
- 26" Ultra-light "Spinergy LX" 02102 **£ 742**
- 24" Off-road (aluminium handrim black) 01392 **£ 164**
- 22" with brake handles for companion (drum brake) (aluminium handrim black) 03385 **£ 284**
- 24" with brake handles for companion (drum brake) (aluminium handrim black) 03384 **£ 284**
- One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°
- ☐ Right ☐ Left 02103 **£ 327**

## Rear wheel tyre

- ☆ High pressure tyre (22"-24"-25"-26")
- ☐ High pressure profile tyre (24")
- Marathon Plus (24"-25") with anti-puncture protection 01669 **£ 55**
- Solid tyre Kik Mako - Colour Black (24") 01041 **£ 87**
- Solid tyre Kik Mako - Colour Grey (24") 01751 **£ 87**

## Handrim - Wheel distance

☆ 3 cm ☐ 4 cm

- ☆ Aluminium (22"-24"-25"-26") 09170 **£ 109**
- Aluminium black SOFT-TOUCH (24"-25") 01684 **£ 109**
- Chromed Steel (24"-25") 00016 **£ 218**
- Titanium (24"-25") 00017 **£ 262**
- Ultra Grepp (24"-25") 00834 **£ 305**
- Max Grepp (24"-25") 00632 **£ 305**
- Ergo Para (24") 06427 **£ 316**
- SURGE LT (24"-25") Oval countour with a rubber insert 00693 **£ 175**
- Quad (24") 04050 **£ 142**
- Silicon cover for standard aluminium handrim Available for measures (22"-24"-25"-26")

## Wheel locks

- ☆ Bent lever push to lock 
- ☐ Straight lever Push to lock 
- ☐ Straight lever Pull to lock
- Sport brake OUT-FRONT  06428 **£ 109**
- Brake with extension lever push to lock  01949 **£ 55**
- Brake with extension lever pull to lock 03372 **£ 55**

English


Date ..... Order ..... User name ..... Address .....

## Footplate - See reference tables N°3

☆ Plastic footrest, rake adjustable

Minimum footplate distance 27cm.


☆2/3 Internal ☐ 2/3 External



☐ One piece footplate - Flip-up! rake and depth adjustable

Minimum footplate distance 26cm.

☐ Internal ☆2/3 Internal ☐ External




○ Minimum footplate distance 21cm. Available only with NEW!!!Multi-function and rotative system - 03945 £ 164,00

☐ One piece aluminium footplate - Flip-up! rake and depth adjustable

Minimum footplate distance 26cm.

☐ Internal ☆2/3 Internal ☐ 2/3 External ☐ External



○ One piece Carbon fibre footplate rake adjustable, swing away, flip-up.

Minimum footplate distance 26cm.

☐ Internal ☐ External



03291 £ 240

## Mudguard


### Mudguard (Aluminium standard frame, Titanium frame TMAX)

☆ ABS mudguard

☐ ABS Straight sideguard

○ Carbon fibre mudguard

○ Carbon fibre Straight sideguard



03946 £ 262

08918 £ 262

### Mudguard (Aluminium frame DYNAMIC)

☆ ABS mudguard

○ Carbon fibre mudguard



03946 £ 262

### Removable version (Aluminium standard frame, Aluminium frame DYNAMIC)

○ Removable ABS mudguard

○ Removable ABS fibre straight sideguard

○ Removable carbon fibre mudguard

○ Removable Carbon fibre straight sideguard



08921 £ 109

08921 £ 109

08919 £ 371

08920 £ 371

## Armrests

○ Swing away armrests

(Swing away, removable and height adjustable armrest)

○ Standard Support (Not with frame Dynamic)

○ Steel support with safety locking system



04828 £ 131

00694 £ 218

### Removable armrest with ABS straight sideguard

○ Desk armrests

(Removable and height adjustable)



00026 £ 164

○ Flip-up armrests

Height and depth adjustable - Swing away and removable



01675 £ 164

○ Removable armrest with ABS straight sideguard



08881 £ 196

○ Removable armrest with carbon-fiber straight sideguard

08882 £ 393

### Kind of pad to be mounted on chosen armrest

☆ Pad STANDARD

○ Ergonomic pad (30x6 cm.)

04104 £ 55

English

Date ..... Order ..... User name ..... Address .....

## SAFETY

☐ 4 points tie-down system rings



08654 £ 164

## Optional 1

☐ Transit wheels



00836 £ 131

☐ Anti-tipper with rapid hookup

(inside rotation to reduce the encumbrance)

☐ Right  
☐ Left



00023 £ 109  
00024 £ 109

☐ Tipping aid

☐ Right  
☐ Left



00025 £ 55  
01046 £ 55

☐ Crutches support

☐ Right  
☐ Left



00838 £ 65  
00837 £ 65

## Optional 2

☐ Quick release front forks

00030 £ 131

☐ Tetraclick

(Quick release axle for quadriplegic)

03947 £ 131

☐ Pad for front frame protectors

07058 £ 39

☐ Rear wheels plate extension

(Not with frame Dynamic)

01950 £ 131

☐ Handybike adaptor

00506 £ 55

☐ Push handles for companion

(Removable and height adjustable)



01048 £ 153

☐ Backrest with folding handles



02958 £ 218

☐ Harness

04103 £ 109

☐ Safety pelvic belt

03972 £ 55

☐ Detachable table

00508 £ 240

☐ Lumbar support

(to insert in the backrest upholstery)

00013 £ 33

☐ Super Head rest

(Folding with the wheelchair)



04049 £ 360

☐ Swing away lateral support

Not available with Posture backrest PHYSIO

☐ Right  
☐ Left



02097 £ 218  
02098 £ 218

☐ Spoke guards Black

(22"-24")

01043 £ 131

☐ Spoke guards Transparent

(22"-24")

00835 £ 131

☐ Spoke guards Progeo design

(22"-24")

00019 £ 153

☐ Spoke guards Dragon

(22"-24")

02106 £ 153

## Additional Kit

☐ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378 £ 436

☐ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357 £ 764

English

Date ..... Order ..... User name ..... Address .....

**PHYSIO - Posture backrest**

☐ Aluminum shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41  
Height ☐ 52

05308 £ 884  
05309 £ 938



Lateral retaining wing ☐ A=high ☐ B=low

☐ Aluminum shell - Lateral retaining DEEP 15cm

Height ☐ 25 ☐ 33 ☐ 41

05310 £ 938

**H25 not available for seat width 39-42-45-48**



Height ☐ 52

05311 £ 993

Lateral retaining wing ☐ A=high ☐ B=low

☐ Carbon shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41

05312 £ 1015



Lateral retaining wing ☐ A=high ☐ B=low

**COVER**


☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems



☐ SWING AWAY - Lateral support

**Not available with Physio backrest with lateral retaining DEEP and in carbon backrest**

**Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)**



Large (17,5x11,5)


☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

**Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)**




Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

**Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)**



Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

(with additional kit)

£ 175  
£ 175

£ 175  
£ 175

☐ Headrest adjustable - Double ball mounting system

**Not available for Physio backrest's height 25 and with carbon shell**



☐ Large (34x14)  
☐ Small (28x14)

05648 £ 393

English

Date ..... Order ..... User name ..... Address .....

Rehateam s.r.l. Vico Negrilli 5 - 31040 Castagnole di Paese (TV) - Tel. +39.0422.484657 / +39.0422.480387 - Fax +39.0422.484661 / +39.0422.484684  
P.IVA IT 03182520266 - C.C.I.A.A. TV n.225867 - Iscr.Trib. TV n.41690 Capitale Sociale I.V. Euro 100.000,00

Email: info@rehateamprogeo.com - web: http://www.rehateamprogeo.com

## Note

Italia Certified company normative UNI EN ISO 9001

Important notices page only if any optional have been chosen

Pag. 6 / 6




Wheelchair : EXELLE VARIO - Rev. 2019.0





Prices in EURO - VAT excluded  
to be added to the price list

☆	standard	(supplied if no requests are indicated)
□	on request	(without extra charge)
○	optional	(cost to be added at the basic price)

## Enclosure to wheelchair order form Exelle Vario

### REFERENCE TABLES

1	Front height				
With closed frame (95°) max 5" caster					
FORK	Casters				
	3"	4"	5"	6"	7"
SMALL 	min. 39	min. 40	min. 42		
	max 45	max 46	max 48		
MEDIUM 	min. 45	min. 45	min. 45	min. 46	min. 48
	max 51	max 52	max 53	max 54	max 55
LONG 	min. 48	min. 48	min. 48	min. 48	min. 48
	max 54	max 55	max 56	max 57	max 58
The size of the fork will be defined only by Rehateam Progeo in accordance with wheelchair final configuration.					

2	Rear height				
Rear frame		Rear wheels			
		22"	24"	25"	26"
STANDARD 	Straight plate 	min. 37	min. 39	min. 40	min 41
		max 45	max 47	max 48	max 49
	Shaped plate 	min. 36	min. 38	min. 39	min. 40
		max 49	max 51	max 52	max 53
DYNAMIC 		min. 36	min. 38	min. 39	min. 40
		max 49	max 51	max 52	max 53

3 Footplate option		Footplate distance		
		Frame Angle 110°	Frame Angle 120°	Frame Angle 95°
Plastic footrest		min. 27	min. 27	min. 27
One piece footplate - Flip-up!		min. 24	min. 26	min. 24
One piece footplate - Flip-up! , with multifunction rotative system		min. 19	min. 21	min. 19
One piece aluminium footplate - Flip-up!		min. 24	min. 26	min. 24
One piece Carbon fibre footplate		min. 24	min. 26	min. 24
Calculated with the best configuration and with the footplates positioned on inside.				
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)				

4 Rear wheels	
Model	Features
22" Ultra-light Progeo - Standard!	(not available in anthracite color)
24" Ultra-light Progeo - Standard!	(available for all couloors)
24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
26" Progeo - Standard!	(colour: black)
26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)