




# Free Form Seating

custom seating. simplified.

Free Form Seating provides a custom fit seating system that can change with a person's postural support needs. This adaptive modular system gives the fitter the ability to contour the shell to the exact needs of the user. As a person's weight fluctuates, as young people grow and as people develop a need for more or less support, this custom seating solution will adapt.

Free Form Seating can be grown, reused, recycled and continuously adapted. Custom seating, simplified.





## Custom back in a box

three back kit sizes. infinite possibilities.

In the box:

- Free Form Shell
- Cover (Air Mesh and foam)
- Mounting hardware
- Aluminum Stays
- 5 Spare Free Form modules
- Extended screws for mounting accessories
- Hex keys

## Why it works:

- Allows for gradual corrective positioning
- Adapts to weight fluctuations
- Grows with children (and post surgical)
- Lightweight, breathable and strong
- Adapters for many types of wheelchairs and strollers
- Heat and moisture dissipation directly through the shell
- Modules can be removed or recessed for total off-loading
- Free Form Shell can be covered in the most appropriate material for the client
- Connecting units (2-balls & edge balls) dampen vibration and give some flex for dynamic movement
- Available in Back Kits, custom shapes or sheets for assembly into any shape

