

Date .....

User name .....

Company stamp

Responsible .....

## Wheelchair measures

Seat width (LS) ☐33 ☐36 ☐39 ☐42 ☐45 ☐48  
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.

Seat depth (PS) ☐35 ☐37,5 ☐40 ☐42,5 ☐45 ☐47,5 ☐50  
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.


Footplate distance (DP) ..... (Real distance between seat and footplate)  
**IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)**

Setting ☒ Standard ☐ Prudential ☐ Active ☐ Extreme

Backrest Height (TS) Std.Value 33,5 .....(min 24.5-max 47)

Camber ☒ 0° ☐ 2° ☐ 4°

### ★ Standard rear adjustment



- ★ Carbon fibre mudguard
- ☐ Carbon fibre Straight sideguard
- ☐ Removable version




08921

£ 109

☐ FULL CARBON rear adjustment plate with integrated mudguard  
(Lightest configuration)

07646

£ 818



- ★ Version mudguard (AVAILABLE ONLY WITH REAR WHEELS 24")
- ☐ Version straight sideguard

★ Polished transparent painting



☐ Opaque transparent painting

08923

£ 164

### Front height - Tick the box with the required height

The possible difference in height between HIGH frame and LOW frame is 3cm (with same fork, caster and assembly hole).

	 SPORT fork (H95)	 Small size fork (H115)
<b>LOW FRAME 90°</b> Front height (min 44,5 - max 50) Footplate distance (min 33)  Footplate distance (min 21) with U-TUBE Aluminium footplate	<p>Casters</p> <p>3" <input type="checkbox"/> ~44,5 <input type="checkbox"/> ~45,5</p> <p>4" <input type="checkbox"/> ~45,5 <input type="checkbox"/> ~46,5</p> <p>5" <input type="checkbox"/> ~48</p>	<p>Casters</p> <p>3" <input type="checkbox"/> ~46,5 <input type="checkbox"/> ~47,5</p> <p>4" <input type="checkbox"/> ~46,5 <input type="checkbox"/> ~47,5 <input type="checkbox"/> ~48,5</p> <p>5" <input type="checkbox"/> ~49 <input type="checkbox"/> ~50</p>
<b>HIGH FRAME 90°</b> Front height (min 47,5 - max 53) Footplate distance (min 36)  Footplate distance (min 24) with U-TUBE Aluminium footplate	<p>Casters</p> <p>3" <input type="checkbox"/> ~47,5 <input type="checkbox"/> ~48,5</p> <p>4" <input type="checkbox"/> ~48,5 <input type="checkbox"/> ~49,5</p> <p>5" <input type="checkbox"/> ~51</p>	<p>Casters</p> <p>3" <input type="checkbox"/> ~49,5 <input type="checkbox"/> ~50,5</p> <p>4" <input type="checkbox"/> ~49,5 <input type="checkbox"/> ~50,5 <input type="checkbox"/> ~51,5</p> <p>5" <input type="checkbox"/> ~52 <input type="checkbox"/> ~53</p>

The real front-seat height could be slightly different from the chosen one considering it is affected by seat depth and rear seat height.

### Rear height .....

The difference between Font Height and Rear Height must be in the 3 - 12 cm range.

- Rear wheels 22" (min 36 - max 44)
- Rear wheels 24" (min 38 - max 46)
- Rear wheels 25" (min 39 - max 47)
- Rear wheels 26" (min 41 - max 48)


Date ..... Order ..... English ..... User name ..... Address .....

**Seat**

- ☆ Seat cloth with pocket (under seat)
- ☐ Seat with tension straps 08883 **£ 196**
- ☐ Carbon rigid seat 08336 **£ 415**
- ☐ Pocket (under seat) [25x11x4,5] 08281 **£ 65**

**CUSTOMIZE COLORS**

- 1 Frame colour**
- ☐ Pearl white ☐ Opaque orange
- ☐ Opaque black 08341 **£ 55**
- ☐ Bright silver
- ☐ Bright red
- ☐ Bright blue
- 2 Fork colour**
- ☆ Opaque black
- ☐ Bright silver ☐ Bright blue
- ☐ Pearl white ☐ Bright fuchsia
- ☐ Bright red ☐ Bright lime
- 08342 **£ 55**
- 3 Colour rear wheels**
- To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)
- ☆ Black
- ☐ Bright silver ☐ Bright blue
- ☐ Pearl white ☐ Bright fuchsia
- ☐ Bright red ☐ Bright lime
- 04692 **£ 55**



**Rear wheels axle**

- ☆ Aluminum rear wheels axle
- ☐ Carbon fiber rear wheels axle (Lightest configuration) **07647 £ 131**


**Backrest with rake adjustable from -4° to +16°**

- ☆ Adjustable velcro straps and padded transpiring upholstery
- ☐ Light adjustable velcro straps and padded transpiring upholstery SPORT 08337 **£ 65**
- ☆ Backrest with double locking (upright and folded)
- ☐ Backrest with upright locking system ONLY (no locking when folded)
- ☐ Fixed back(non-folding back)
- ☆ Aluminium backrest
- ☐ Titanium backrest
- ☐ Carbon fibre backrest (Lightest configuration) **06811 £ 109**
- ☐ A I R - Integrated posture backrest (not removable) - Lateral retaining 6,5cm. 04389 **£ 600**
- not available for seat width 48 - Std. Backrest upholstery not provided - Handle tube aluminium with short tubes 06815 **£ 1015**



**Lateral retaining wing B=low**

**Carbon shell H25**



Backrest Height (TS) min 29 ..... max 41


**Carbon shell H33**



Backrest Height (TS) ..... (37 + 46)

**Lateral retaining wing A=high**

**Carbon shell H30**



Backrest Height (TS) ..... (32 + 42,5)


**AIR Cover**

- ☆ Cover Standard - fabric insert: GREY, stitchings GREY

- ☐ **Cover SPORT** 09195 **£ 87**
- ☐ fabric insert: BLACK - Stitchings GREY
- ☐ fabric insert: BLACK - Stitchings LIME
- ☐ fabric insert: BLACK - Stitchings ORANGE
- ☐ fabric insert: BLACK - Stitchings RED
- ☐ fabric insert: BLACK - Stitchings BLUE
- ☐ fabric insert: BLACK - Stitchings BLACK
- ☐ fabric insert: BLACK - Stitchings FUCHSIA
- ☐ fabric insert: BLACK - Stitchings ANTHRACITE

- ☐ **Cover ACTIVE** 09195 **£ 87**
- ☐ fabric insert: RED - Stitchings BLACK
- ☐ fabric insert: LIME - Stitchings BLACK
- ☐ fabric insert: BLUE - Stitchings BLACK

fabric insert:



Transpiring fabric stitchings

English

Date Order User name Address

### Handle tube

☆ Handle tube STD

- ☐ Handle tube SPORT Without push handles  
☐ Handle tube with short tubes

☐ Carbon fibre handle (with short tubes) Weight reduction approx. 110gr.  
**Not available with: Posture backrest PHYSIO - Posture backrest A I R**

04390 £ 109

☐ Folding handles aluminium



02958 £ 218

### Casters

☆ 4" Soft

- ☐ 3" Sport  
☐ 5" Soft

- ☐ 3" Active Soft (wide) (Aluminium alloy rim)  
☐ 4" Active Sport (narrow) (Aluminium alloy rim)  
☐ 4" Active Soft (wide) (Aluminium alloy rim)  
☐ 5" Active Soft (wide) (Aluminium alloy rim)

06824 £ 109

01398 £ 109

03943 £ 109

07518 £ 109

### Rear wheels

☆ 24" Ultra-light Progeo - Standard!

- ☐ 22" Ultra-light Progeo - Standard!  
☐ 25" Ultra-light Progeo - Standard!  
☐ 26" Progeo - Standard!

(colour: black)

(colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT

Weight reduction approx. 450gr.

04383 £ 742

☐ 24" Ultra-light "Spinergy LX"

01685 £ 742

☐ 24" Ultra-light "Spinergy LXL"

☐ Black hub and white blade spokes

☐ Hub and blade spokes black

07846 £ 742

☐ 24" Ultra-light "Spinergy Flex rim"

02270 £ 1200

☐ 25" Ultra-light "Spinergy LX"

02101 £ 742

☐ 25" Ultra-light "Spinergy LXL"

☐ Black hub and white blade spokes

07847 £ 742

☐ 26" Ultra-light "Spinergy LX"

02102 £ 742

☐ 24" Off-road (aluminium handrim black)

01392 £ 164

### Rear wheel tyre

☆ High pressure tyre (22"-24"-25"-26")

☐ High pressure profile tyre (24")

☐ Marathon Plus (24"-25")

with anti-puncture protection

01669 £ 55

☐ Solid tyre Kik Mako - Colour Black (24")

01041 £ 87

☐ Solid tyre Kik Mako - Colour Grey (24")

01751 £ 87

### Handrim - Wheel distance

☆ 3 cm □ 4 cm

☆ Aluminium (22"-24"-25"-26")

09170 £ 109

☐ Aluminium black SOFT-TOUCH (24"-25")

01684 £ 109

☐ Chromed Steel (24"-25")

00016 £ 218

☐ Titanium (24"-25")

00017 £ 262

☐ Ultra Grepp (24"-25")

00834 £ 305

☐ Max Grepp (24"-25")

00632 £ 305

☐ Ergo Para (24")

06427 £ 316

☐ SURGE LT (24"-25")

Oval countour with a rubber insert

00693 £ 175

☐ Quad (24") (24")

04050 £ 142

☐ Silicon cover for standard aluminium handrim Available for measures (22"-24"-25"-26")

### Wheel locks

☆ Bent lever push to lock

- ☐ Straight lever Push to lock  
☐ Straight lever Pull to lock

☐ Sport brake OUT-FRONT

06428 £ 109

☐ Brake with extension lever push to lock

01949 £ 55

☐ Brake with extension lever pull to lock

03372 £ 55



Date ..... Order ..... English ..... User name ..... Address .....

## Footplate

☆ Aluminium footplate. Rake not adjustable

☆ Aluminium footplate ☐ with plastic plate

☐ One piece aluminium footplate - Flip-up! FLIP-BACK! - rake and depth adjustable

☐ Aluminium footplate - rake and depth adjustable



☆ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External

☐ U-TUBE - Aluminium footplate, rake and depth adjustable

Low Frame minimum footplate distance 21 - High Frame minimum footplate distance 24

☐ One piece Carbon fibre footplate. Rake adjustable

☆ Internal ☐ External




09420 £ 164



04456 £ 218



09853 £ 196



04110 £ 327

## Armrests

☐ Swing away armrests (Swing away, removable and height adjustable armrest)



00694 £ 218

☐ Removable armrest with carbon-fiber straight sideguard  
Available only with: Standard rear adjustment



09079 £ 196

### Kind of pad to be mounted on chosen armrest

☒ Pad STANDARD  
☐ Ergonomic pad (30x6 cm.)

04104 £ 55

## Optional 1

### Only one optional devices for side can be chosen

☐ Transit wheels Exclude the installation of Anti-tippers



03290 £ 131

☐ Anti-tipper with rapid hookup  
 Exclude the installation of Transit wheels  
 (inside rotation to reduce the encumbrance)  
☐ Right  
☐ Left



06821 £ 109

06822 £ 109

## Optional 2

☐ Push handles for companion (Removable and height adjustable)



03645 £ 153

☐ Pad for front frame protectors

07058 £ 39

☐ Tetraclip (Quick release axle for quadriplegic)

03947 £ 131

☐ Safety pelvic belt

03972 £ 55

☐ Lumbar support (to insert in the backrest upholstery)

00013 £ 33

☐ Swing away lateral support  
 Not available with:  
 Posture backrest PHYSIO - Posture backrest A I R

☐ Right  
☐ Left



02097 £ 218

02098 £ 218

☐ Spoke guards Black (22"-24")  
☐ Spoke guards Transparent (22"-24")  
☐ Spoke guards Progeo design (22"-24")  
☐ Spoke guards Dragon (22"-24")

01043 £ 131

00835 £ 131

00019 £ 153

02106 £ 153

## Additional Kit

☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour)

03378 £ 436

☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour)

09357 £ 764

Date ..... Order ..... English User name ..... Address .....

**PHYSIO - Posture backrest**

- ☐ Aluminum shell - Lateral retaining STANDARD 7,5cm Height ☐ 25 ☐ 33 ☐ 41 05308 £ 884  
Height ☐ 52 05309 £ 938



Lateral retaining wing ☐ A=high ☐ B=low

- ☐ Aluminum shell - Lateral retaining DEEP 15cm Height ☐ 25 ☐ 33 ☐ 41 05310 £ 938  
**H25 not available for seat width 39-42-45-48**



Height ☐ 52

05311 £ 993

Lateral retaining wing ☐ A=high ☐ B=low

- ☐ Carbon shell - Lateral retaining STANDARD 7,5cm Height ☐ 25 ☐ 33 ☐ 41 05312 £ 1015



Lateral retaining wing ☐ A=high ☐ B=low


**COVER**

☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




- ☐ SWING AWAY - Lateral support  
Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width,  
adjustable (0 ÷ -1,5cm)




- Large (17,5x11,5)  
☐ Left 05646 ☐ Right 05647  
Small (14x9,5)  
☐ Left 05642 ☐ Right 05643  
☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

Internal thoracic support to the backrest width,  
adjustable (-4 ÷ -5,5cm)



- Large (17,5x11,5)  
☐ Left 05646 ☐ Right 05647  
Small (14x9,5)  
☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width,  
adjustable (-8,5 ÷ -10cm)



- Large (17,5x11,5)  
☐ Left 05646 ☐ Right 05647 £ 175  
Small (14x9,5)  
☐ Left 05642 ☐ Right 05643 £ 175  
(with additional kit)

- ☐ Headrest adjustable - Double ball mounting system  
Not available for Physio backrest's height 25 and with carbon shell



- ☐ Large (34x14) 05648 £ 393  
☐ Small (28x14)

English

Date ..... Order ..... User name ..... Address .....

Rehateam s.r.l. Vicolo Negrilli, 5 - 31040 Castagnole di Paese (TV) - Tel. +39.0422.484657 / +39.0422.480387 - Fax +39.0422.484661 / +39.0422.484684  
P.IVA IT 03182520266 - C.C.I.A.A. TV n.225867 - Iscr. Trib. TV n.41690 Capitale Sociale I.V. Euro 100.000,00

Email: info@rehateamprogeo.com - web: http://www.rehateamprogeo.com

1	<b>Rear wheels</b>	
	<b>Model</b>	<b>Features</b>
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
	26" Progeo - Standard!	(colour: black)
26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)	
<b>Note</b>		
Italia Certified company normative UNI EN ISO 9001		Important notice! s page only if any optional have been chosen
Pag. 6 / 6	Wheelchair : JOKER R2 - Rev. 2019.0	Prices in EURO - VAT excluded to be added to the price list
		<input checked="" type="checkbox"/> standard (supplied if no requests are indicated) <input type="checkbox"/> on request (without extra charge) <input type="checkbox"/> optional (cost to be added at the basic price)