


Date

User name

Company stamp

Responsible

Wheelchair measures							
Seat width (LS)	<input type="checkbox"/> 33	<input type="checkbox"/> 36	<input type="checkbox"/> 39	<input type="checkbox"/> 42	<input type="checkbox"/> 45	<input type="checkbox"/> 48	
Seat depth (PS)	<input type="checkbox"/> 35	<input type="checkbox"/> 37,5	<input type="checkbox"/> 40	<input type="checkbox"/> 42,5	<input type="checkbox"/> 45	<input type="checkbox"/> 47,5	<input type="checkbox"/> 50
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.							
Footplate distance (DP)	(Real distance between seat and footplate)						
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)							
Setting	☆ Standard		<input type="checkbox"/> Prudential	<input type="checkbox"/> Active	<input type="checkbox"/> Extreme		
Backrest Height (TS)	Std.Value 33,5(min 24.5-max 47)						
camber	☆ 0° <input type="checkbox"/> 2° <input type="checkbox"/> 4°						
Front height STD 49			Rear height STD 42		



☆ Alluminium rear adjustment plate

☆ Removable mudguard ABS

☐ Removable ABS fibre straight sideguard

☐ Removable carbon fibre mudguard


☐ Removable Carbon fibre straight sideguard

03946

£ 262

08918

£ 262



☐ FULL CARBON rear adjustment plate with integrated mudguard

07646

£ 818

☆ Version mudguard (Available only with rear wheels 24")

☐ Version straight sideguard

☆ Polished transparent painting

☐ Opaque transparent painting

08923

£ 164

☐ CARBON PACK


(Lightest configuration)

07650

£ 1909

FULL CARBON rear adjustment plate with integrated mudguard - Carbon fibre frame - Carbon fiber rear wheels axle

Savings £ 349 on same optional total extracharge



☆ Version mudguard (AVAILABLE ONLY WITH REAR WHEELS 24")

☐ Version straight sideguard


☆ Polished transparent painting

☐ Opaque transparent painting

07752

£ 273

JOKER



☆ Aluminium footplate, rake adjustable

Minimum footplate distance 35cm.

☆ Internal ☐ 2/3 Internal

☐ External ☐ 2/3 External

☐ Aluminium footplate, with titanium tubes - rake and depth adjustable

Minimum footplate distance 35cm.

☆ Internal ☐ 2/3 Internal

☐ External ☐ 2/3 External

☐ One piece Carbon fibre footplate


Minimum footplate distance 35cm.

not available for seat width 48

☐ One piece Carbon fibre footplate, rake adjustable

Minimum footplate distance 35cm.

☆ Internal ☐ External




07649

£ 87



03348


£ 131



03346

£ 218

JOKER ZERO



☐ Plastic footrest, rake adjustable

Minimum footplate distance 15cm.

☆ 2/3 Internal ☐ 2/3 External

☐ One piece aluminium footplate - Flip-up!, rake and depth adjustable

Minimum footplate distance 15cm.

☐ Internal ☐ 2/3 Internal

☐ External ☐ 2/3 External

☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up.

Minimum footplate distance 15cm.

☆ Internal ☐ External




06814

£ 164



06856

£ 164



06857

£ 349

Date	Order	English User name	Address
------------	-------------	----------------------------	---------------

Handle tube			
☆ Handle tube STD			
<input type="checkbox"/> Handle tube with short tubes			
<input type="checkbox"/> Handle tube SPORT Without push handles			
<input type="checkbox"/> Carbon fibre handle (with short tubes)	(Lightest configuration)	04390	£ 109
Not available with: Posture backrest PHYSIO - Posture backrest A I R			
<input type="checkbox"/> Folding handles aluminium		02958	£ 218



Casters			
☆ 4" Soft			
<input type="checkbox"/> 3" Sport			
<input type="checkbox"/> 5" Soft			
<input type="checkbox"/> 6" Soft			
<input type="checkbox"/> 6" Pneumatic			
(Lightest configuration)			
<input type="checkbox"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824	£ 109
<input type="checkbox"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	£ 109
<input type="checkbox"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943	£ 109
<input type="checkbox"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518	£ 109

Support front fork			
☆ Front fork support V Design			
<input type="checkbox"/> Front fork support with integrated spirit level			



Special Options			
<input type="checkbox"/> Front fork FROG LEGS	(colour: black)	06426	£ 436
To be matched only with: casters 3" 4" 5"			



Colour set: Rear wheels and forks - To be matched only with:Rear wheels Ultra-light Progeo - Standard! (22 - 24)			
☆ Black 22" 24" 25" 26"			
<input type="checkbox"/> Bright silver	22" 24"	00833	£ 109
<input type="checkbox"/> Bright anthracite	24"	00833	£ 109
<input type="checkbox"/> Bright orange	22" 24"	00833	£ 109
<input type="checkbox"/> Bright fuchsia	22" 24"	00833	£ 109
<input type="checkbox"/> Bright blue	22" 24"	00833	£ 109
<input type="checkbox"/> Bright red	22" 24"	00833	£ 109
<input type="checkbox"/> Bright lime	22" 24"	00833	£ 109
<input type="checkbox"/> Pearl white	24"	00833	£ 109

Rear wheels			
☆ 24" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 26" Progeo - Standard!			
<input type="checkbox"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	Weight reduction approx. 450gr.	04383	£ 742
<input type="checkbox"/> 24" Ultra-light "Spinergy LX"		01685	£ 742
<input type="checkbox"/> 24" Ultra-light "Spinergy Flex rim"		02270	£ 1200
<input type="checkbox"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes <input type="checkbox"/> Hub and blade spokes black	07846	£ 742
<input type="checkbox"/> 25" Ultra-light "Spinergy LX"		02101	£ 742
<input type="checkbox"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes	07847	£ 742
<input type="checkbox"/> 26" Ultra-light "Spinergy LX"		02102	£ 742
<input type="checkbox"/> 24" Off-road (aluminium handrim black)		01392	£ 164

Rear wheel tyre			
☆ High pressure tyre	(22"-24"-25"-26")	(Lightest configuration)	
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="checkbox"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669 £ 55
<input type="checkbox"/> Solid tyre Kik Mako - Colour Black	(24")		01041 £ 87
<input type="checkbox"/> Solid tyre Kik Mako - Colour Grey	(24")		01751 £ 87

Handrim - Wheel distance			
☆ 3 cm		□ 4 cm	
☆ Aluminium	(22"-24"-25"-26")	(Lightest configuration)	
<input type="checkbox"/> Aluminium black SOFT-TOUCH	(24"-25")		09170 £ 109
<input type="checkbox"/> Chromed Steel	(24"-25")		01684 £ 109
<input type="checkbox"/> Titanium	(24"-25")		00016 £ 218
<input type="checkbox"/> Ultra Grepp	(24"-25")		00017 £ 262
<input type="checkbox"/> Max Grepp	(24"-25")		00834 £ 305
<input type="checkbox"/> Ergo Para	(24")		00632 £ 305
<input type="checkbox"/> SURGE LT	(24"-25")	Oval countour with a rubber insert	06427 £ 316
<input type="checkbox"/> Quad (24")	(24")		00693 £ 175
<input type="checkbox"/> Silicon cover for standard aluminium handrim		Available for measures (22"-24"-25"-26")	04050 £ 142

Date Order English User name Address

Wheel locks

☆ Bent lever push to lock

- ☐ Straight lever Push to lock
☐ Straight lever Pull to lock

☐ Sport brake OUT-FRONT

(Lightest configuration)

- ☐ Brake with extension lever push to lock
☐ Brake with extension lever pull to lock



06428

£ 109

01949

£ 55

03372

£ 55

Armrests

☐ Swing away armrests

(Swing away, removable and height adjustable armrest)

☆ Pad STANDARD

☐ Ergonomic pad (30x6 cm.)



00694

£ 218

04104

£ 55

SAFETY

☐ 4 points tie-down system rings



08654

£ 164

Optional 1

Only one optional devices for side can be chosen

☐ Transit wheels

Exclude the installation of Anti-tippers



03290

£ 131

☐ Anti-tipper BASE

Exclude the installation of Transit wheels

☐ Right

☐ Left

06819

£ 87

06820

£ 87

☐ Anti-tipper with rapid hookup
 (inside rotation to reduce the encumbrance)

Exclude the installation of Transit wheels

☐ Right

☐ Left



06821

£ 109

06822

£ 109

Optional 2

☐ Push handles for companion

(Removable and height adjustable)



03645

£ 153

☐ Tetraclip

(Quick release axle for quadriplegic)

03947

£ 131

☐ Pad for front frame protectors

07058

£ 39

☐ Safety pelvic belt

03972

£ 55

☐ Lumbar support

(to insert in the backrest upholstery)

00013

£ 33

☐ Swing away lateral support

Not available with: Posture backrest PHYSIO - Posture backrest A I R

☐ Right

☐ Left



02097

£ 218

02098

£ 218

☐ Spoke guards Black

(22"-24")

01043

£ 131

☐ Spoke guards Transparent

(22"-24")

00835

£ 131

☐ Spoke guards Progeo design

(22"-24")

00019

£ 153

☐ Spoke guards Dragon

(22"-24")

02106

£ 153

Additional Kit

☐ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378

£ 436

☐ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357

£ 764

Date Order English User name Address

PHYSIO - Posture backrest

☐ Aluminum shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41
Height ☐ 52

05308 £ 884
05309 £ 938



Lateral retaining wing ☐ A=high ☐ B=low

☐ Aluminum shell - Lateral retaining DEEP 15cm

Height ☐ 25 ☐ 33 ☐ 41

05310 £ 938

H25 not available for seat width 39-42-45-48



Height ☐ 52

05311 £ 993

Lateral retaining wing ☐ A=high ☐ B=low

☐ Carbon shell - Lateral retaining STANDARD 7,5cm

Height ☐ 25 ☐ 33 ☐ 41

05312 £ 1015



Lateral retaining wing ☐ A=high ☐ B=low

COVER


☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems



☐ SWING AWAY - Lateral support

Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)



Large (17,5x11,5)


☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)




Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)



Large (17,5x11,5)

☐ Left 05646 ☐ Right 05647

Small (14x9,5)

☐ Left 05642 ☐ Right 05643

(with additional kit)

£ 175
£ 175

£ 175
£ 175

☐ Headrest adjustable - Double ball mounting system

Not available for Physio backrest's height 25 and with carbon shell



☐ Large (34x14)
☐ Small (28x14)

05648 £ 393

English

Date Order User name Address

Rehateam s.r.l. Vicolo Negrelli, 5 - 31040 Castagnole di Paese (TV) - Tel. +39.0422.484657 / +39.0422.480387 - Fax +39.0422.484661 / +39.0422.484684
P.IVA IT 03182520266 - C.C.I.A.A. TV n.225867 - Iscr.Trib. TV n.41690 Capitale Sociale I.V. Euro 100.000,00

Email: info@rehateamprogeo.com - web: http://www.rehateamprogeo.com

Note

Italia Certified company normative UNI EN ISO 9001

Important notice! s page only if any optional have been chosen

Pag. 6 / 6




Wheelchair : JOKER - Rev. 2019.0

Prices in EURO - VAT excluded
to be added to the price list

☒ standard (supplied if no requests are indicated)
☐ on request (without extra charge)
☐ optional (cost to be added at the basic price)

Enclosure to wheelchair order form Joker

REFERENCE TABLES

1	Front height				
With closed frame (95°) max 5" caster					
With closed frame (90°) max 5" caster					
With closed frame (85°) max 4" caster					
FORK		Casters			
		3"	4"	5"	6"
SPORT		max 52	max 53	max 54	
SMALL		max 54	max 55	max 56	max 57
MEDIUM		max 56	max 57	max 58	max 59
The size of the fork will be defined only by Rehateam Progeo in accordance with wheelchair final configuration.					

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min	min	min	min	
36	38	39	40	
max	max	max	max	
44	46	47	48	

3	Footplate option	Footplate distance
	Aluminium footplate rake adjustable	min 35
	One piece Carbon fibre footplate	min 35
	One piece Carbon fibre footplate rake adjustable	min 35
	Plastic footrest, rake adjustable	min 35
	One piece aluminium footplate - Flip-up!, rake and depth adjustable	min 35
	One piece Carbon fibre footplate rake adjustable, swing away, flip-up.	min 35
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)		
Calculated with the best configuration and with the footplates positioned on inside.		
NOTE! Lower footplate distances are possible only after our valuation		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)