

The compact leg and upper body exerciser with motor

# mobi 540



**THERA-soft®**  
connectable + 

● = Equipment variant ○ = Accessory

THERA-Trainer mobi 540	
Item no.	A004-419
<b>NZD excl Freight and GST</b>	<b>\$ 6495</b>

Name	Item no.	Equipment	COST
Base unit	A003-975	●	
Foot rests	A003-974	●	
Therapy grips	A001-420	●	
Control and display unit mobi 2.7" wireless	A006-172	●	
Wristbands for therapy grip/safety handle	A002-440	○	POA
Carry bag	A003-976	○	POA
THERA-soft	A006-630	○	POA
Media station	A008-119	○	POA

The compact leg and upper body exerciser with motor

# mobi 540

Product specifications	
<b>Product description</b>	THERA-Trainer mobi is a powered training device using cyclical movement actively and/or passively mobilizing the upper or lower extremities to redevelop muscles, restore motion to joints and improve the cardiovascular function.
<b>Special product features</b>	<ul style="list-style-type: none"> <li>- <b>Foot rests with safety edge and velcro straps to secure the foot:</b> Ensure a safe workout through optimal support and without foot slippage.</li> <li>- <b>Very good stability and slip resistance:</b> Enables safe training without constant slipping.</li> <li>- <b>Quick and easy cleaning:</b> Thorough wipe disinfection is possible thanks to smooth surfaces.</li> <li>- <b>Tool-free replacement of foot rests and therapy grips:</b> Enables quick and easy adaptation to different treatment settings and therapy goals.</li> <li>- <b>Forward-looking treatment:</b> The THERA-Trainer mobi can be integrated into a THERA-Trainer complete solution at any time.</li> </ul>
<b>Place of use</b>	Suitable for use in enclosed spaces in facilities e.g. clinics, nursing homes, therapy practices or at the patient's home
<b>Technical specification</b>	<ul style="list-style-type: none"> <li>- <b>Dimensions (LxWxH):</b> approx. 50 cm x 46 cm - 56 cm x 43 cm</li> <li>- <b>Weight:</b> approx. 13.2 kg (with foot rests) or 12.1 kg (with therapy grips)</li> <li>- <b>Electric motor:</b> Nominal power = 100 watts, continuous power = 50 watts</li> <li>- <b>Power connection:</b> 1x 100-240 V, 50/60 Hz</li> <li>- <b>Protection class:</b> II</li> <li>- <b>Degree of protection:</b> Type BF</li> <li>- <b>Protection category:</b> IP 21</li> <li>- <b>Ambient conditions when in use:</b>            Temperature: 5 °C to 40 °C            Relative humidity: 15 % to 93 % Rh            Air pressure: 700 hPa to 1060 hPa</li> <li>- <b>Sound emission:</b> LpA &lt; 70 dB (A)</li> <li>- <b>Service life:</b> 7 years</li> </ul>
<b>Compatibility</b>	THERA-soft
<b>Patient requirements</b>	<ul style="list-style-type: none"> <li>- <b>Patient weight:</b> max. 140 kg</li> <li>- <b>Patient height:</b> 120 - 200 cm</li> </ul>
<b>Delivery and commissioning</b>	<ul style="list-style-type: none"> <li>- <b>Packaging (LxWxH):</b> 56 cm x 46 cm x 40 cm</li> <li>- <b>Delivery:</b> Package via parcel service</li> <li>- <b>Transport and storage conditions:</b>            Temperature: -25 °C to 70 °C            Relative humidity: 15 % to 93 % Rh            Air pressure: 700 hPa to 1060 hPa</li> <li>- <b>Commissioning:</b> by authorised specialists (recommended)</li> </ul>
<b>Service / Maintenance</b>	<ul style="list-style-type: none"> <li>- <b>Service:</b> only by authorised service technicians</li> <li>- <b>Maintenance:</b> maintenance-free</li> <li>- <b>Technical inspection:</b> safety inspection every 24 months (manufacturer recommendation)</li> </ul>
<b>Registration</b>	Medical device class: IIa CE marking: CE0297 (monitored by certification body)