

Date

User name

Responsible


Carbon fibre frame made-to-measure.

Keep attention to all measures.

NOIR is a wheelchair without possibilities of adjustments after its manufacturing.


See reference tables

Take care to write all the requested measures (in centimetres)



Front height (min 44-max 54)

Rear height (min 37-max 48)
calculated with rear wheels 24"



ESSENCE

☆ CARBON Essence
Polished transparent painting

○ OPAQUE Essence
Opaque transparent painting

○ LIME Essence

○ RED Essence

○ BLUE Essence

○ ANTHRACITE + OPAQUE Essence

08925 £ 545

08926 £ 764

08927 £ 764

08928 £ 764

08929 £ 982

Fork colour

☆ Opaque black


☐ Pearl white

Footplate option

☆ Carbon fibre frame with integrated footplate, not adjustable
(Lightest configuration)

☐ Aluminium footplate adjustable distance +/- 1,5cm (rake and depth adjustable)
Minimum footplate distance 35cm.

☆ Internal ☐ 2/3 Internal
☐ External ☐ 2/3 External




English
Date Order User name Address

Seat

- ☆ Seat cloth with pocket (under seat)
☐ Seat with tension straps 08883 £ 196
☐ Carbon rigid seat 08336 £ 415
☐ Pocket (under seat) [25x11x4,5] 08281 £ 65

Backrest with rake adjustable from -4° to +16°


- ☆ Backrest with double locking (upright and folded)
☐ Backrest with upright locking system ONLY (no locking when folded)
☐ Fixed back (non-folding back)
 ☆ Carbon fibre backrest (Lightest configuration)
☐ Titanium backrest



- ☐ **AIR - Integrated posture backrest (not removable) - Lateral retaining 6,5cm.** 06815 £ 1015
 not available for seat width 48 - Std. Backrest upholstery not provided - Handle tube aluminium with short tubes


Lateral retaining wing B=low

Carbon shell H25



Backrest Height (TS) min 29 max 41


Carbon shell H33



Backrest Height (TS) (37 ÷ 46)

Lateral retaining wing A=high

Carbon shell H30



Backrest Height (TS) (32 ÷ 42,5)

AIR Cover

- ☆ Cover Standard - fabric insert: GREY, stitchings GREY

☐ **Cover SPORT**


09195 £ 87

- ☐ fabric insert: BLACK - Stitchings GREY
☐ fabric insert: BLACK - Stitchings LIME
☐ fabric insert: BLACK - Stitchings ORANGE
☐ fabric insert: BLACK - Stitchings RED
☐ fabric insert: BLACK - Stitchings BLUE
☐ fabric insert: BLACK - Stitchings BLACK
☐ fabric insert: BLACK - Stitchings FUCHSIA
☐ fabric insert: BLACK - Stitchings ANTHRACITE

☐ **Cover ACTIVE**

09195 £ 87

- ☐ fabric insert: RED - Stitchings BLACK
☐ fabric insert: LIME - Stitchings BLACK
☐ fabric insert: BLUE - Stitchings BLACK



Handle tube

- ☆ Handle tube with short tubes
☐ Handle tube SPORT Without push handles (Lightest configuration)
☐ Folding handles aluminium



02958 £ 218

Castors

- ☆ 3" Sport (Lightest configuration)
☐ 4" Soft
☐ 3" Active Soft (wide) (Aluminium alloy rim)
☐ 4" Active Sport (narrow) (Aluminium alloy rim)
☐ 4" Active Soft (wide) (Aluminium alloy rim)
☐ 5" Active Soft (wide) (Aluminium alloy rim)

Rear wheels - See chart N°3

- ☆ 24" ULTRA-LIGHT PROGEO KEY SPORT (Lightest configuration)
☐ 24" Ultra-light "Spinergy LX"
☐ 24" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes ☐ Hub and blade spokes black
☐ 25" Ultra-light "Spinergy LX"
☐ 25" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes
☐ 26" Ultra-light "Spinergy LX"
☐ 24" Ultra-light "Spinergy Flex rim" 04353 £ 458

English

Date Order User name Address

Rear wheel tyre

- ☆ High pressure tyre (22"-24"-25"-26") (Lightest configuration)
- ☐ High pressure profile tyre (24")
- ☐ Marathon Plus (24"-25") with anti-puncture protection 01669 £ 55
- ☐ Solid tyre Kik Mako - Colour Black (24") 01041 £ 87
- ☐ Solid tyre Kik Mako - Colour Grey (24") 01751 £ 87

Handrim - Wheel distance

☆ 3 cm

☐ 4 cm

- ☆ Aluminium (22"-24"-25"-26") (Lightest configuration)
- ☐ Aluminium black SOFT-TOUCH (24"-25") 09170 £ 109
- ☐ Chromed Steel (24"-25") 01684 £ 109
- ☐ Titanium (24"-25") 00016 £ 218
- ☐ Ultra Grepp (24"-25") 00017 £ 262
- ☐ Max Grepp (24"-25") 00834 £ 305
- ☐ Ergo Para (24") 00632 £ 305
- ☐ SURGE LT (24"-25") Oval countour with a rubber insert 06427 £ 316
- ☐ Quad (24") (24") 00693 £ 175
- ☐ Silicon cover for standard aluminium handrim Available for measures (22"-24"-25"-26") 04050 £ 142

Wheel locks

- ☆ Without brake
- ☐ Bent lever push to lock
- ☐ Straight lever Push to lock
- ☐ Straight lever Pull to lock
- ☐ Sport brake OUT-FRONT (Lightest configuration)
- ☐ Brake with extension lever push to lock
- ☐ Brake with extension lever pull to lock

Mudguard

- ☆ Carbon fibre mudguard
- ☐ Carbon fibre Straight sideguard ☆ Aligned wheel ☐ Under wheel cm ☐ Above wheel cm
- ☐ Removable carbon fibre mudguard 04109 £ 218
- ☐ Removable Carbon fibre straight sideguard ☆ Aligned wheel ☐ Under wheel cm ☐ Above wheel cm 04109 £ 218

Armrests

- ☐ Swing away armrests (Swing away, removable and height adjustable armrest) 00694 £ 218
- ☆ Pad STANDARD
- ☐ Ergonomic pad (30x6 cm.) 04104 £ 55

Optional 1

- ☐ Pad for front frame protectors 07058 £ 39
- ☐ Transit wheels 03290 £ 131
- ☐ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance)
- ☐ Right 06821 £ 109
- ☐ Left 06822 £ 109

Date Order English User name Address

Optional 2

- ☐ Lumbar support (to insert in the backrest upholstery) 00013 £ 33
- ☐ Tetraclip (Quick release axle for quadriplegic) 03947 £ 131

Additional Kit

- ☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour) 03378 £ 436
- ☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour) 09357 £ 764

PHYSIO - Posture backrest

Std. Backrest upholstery not provided

- ☐ Aluminum shell
Lateral retaining STANDARD 7,5cm
- Height ☐ 25 ☐ 33 ☐ 41 05308 £ 884
Height ☐ 52 05309 £ 938
Lateral retaining wing ☐ A=high ☐ B=low
- ☐ Aluminum shell
Lateral retaining DEEP 15cm
- Height ☐ 25 ☐ 33 ☐ 41 05310 £ 938
H 25 not available for seat width 39-42-45-48
Height ☐ 52 05311 £ 993
Lateral retaining wing ☐ A=high ☐ B=low
- ☐ Carbon shell
Lateral retaining STANDARD 7,5cm
- Height ☐ 25 ☐ 33 ☐ 41 05312 £ 1015
Lateral retaining wing ☐ A=high ☐ B=low


COVER

- ☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




- ☐ SWING AWAY - Lateral support
Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width,
adjustable (0 ÷ -1,5cm)




- Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647
- Small (14x9,5)
☐ Left 05642 ☐ Right 05643
- ☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

Internal thoracic support to the backrest width,
adjustable (-4 ÷ -5,5cm)



- Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647
- Small (14x9,5)
☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width,
adjustable (-8,5 ÷ -10cm)



- Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647 £ 175
- Small (14x9,5)
☐ Left 05642 ☐ Right 05643 £ 175
- (with additional kit) £ 175

- ☐ Headrest adjustable - Double ball mounting system
Not available for Physio backrest's height 25 and with carbon shell

- ☐ Large (34x14)
☐ Small (28x14)



05648 £ 393

Date Order English User name Address

REFERENCE TABLES

IMPORTANT NOTICE Keep attention to all measures.
NOIR is a wheelchair without possibilities of adjustments after its manufacturing.
(1) The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

1	A Footplate width		
DP	35÷39	40÷45	46÷50
LS			
33	21,5	21	20,5
36	24,5	24	23,5
39	27,5	27	26,5
42	30,5	30	29,5
45	33,5	33	32,5
48	36,5	36	35,5

2	Setting	C :
	Extremely prudential	7,5
	Prudential	9
	Standard	10,5
	Active	12
	Extreme	13,5

3	Maximum width
Camber	
0°	LS + 18
2°	LS + 21
4°	LS + 24

Maximum width calculated including 3 cm. distance between tyre and handrim