

Date

Responsible

Distributor's net price

Information on end user

Wheelchair measures

Name

Surname

Phone

Mail

Heightcm

WeightKg

Wheelchair model and brand.....

Seat width (LS)cm

Seat depth (PS)cm

Footplate ☐ Internal ☐ External

Note:

.....

☆ **Standard version (36V battery)**

- 16" wheel
- Tyre suitable for both road and off-road
- Tubular fork (oval section)
- Disk brake (mechanical)
- Display
- Rear gear
- Cruise control
- Motor with "Regen Battery" function
- 50 km autonomy (approximate distance calculated at speed 15km/h)

○ **Version with 48V battery**

(additional characteristics to the Standard Version)

07936 **337,50**

- External control unit
- Speed and motor brake programmable from the display
- 5 gears
- 60 km autonomy (approximate distance calculated at speed 15km/h)

☐ *White*

☐ *Opaque black*

<input type="radio"/> Hydraulic disk brake	07938	135,00
<input type="radio"/> Double caliper disk brake (mechanical)	07939	120,00
<input type="radio"/> Easel	08196	225,00
<input type="radio"/> Safety kit (lights, bell and rear view mirror)	07940	90,00
<input type="radio"/> Carbon fiber mudguard	07942	93,75
<input type="radio"/> Special handlebar (for quadriplegic)	07941	146,25
<input type="radio"/> Basket (with plate and support)	08208	150,00
<input type="radio"/> Additional pair of clamps	07943	142,50
<input type="radio"/> Travel bag	07944	131,25
<input type="radio"/> Waterproof leg cover	07945	22,50
<input type="radio"/> Spare tyre and tube	07946	37,50

Italia Certified company normative UNI EN ISO 9001

Important notice! is page only if any optional have been chosen