

Supplementary Brochure

Providing equipment for greater
independence, dignity and choice



Neater Solutions for everyday needs

Neater Solutions is renowned worldwide for the design, manufacture and supply of equipment that gives people greater independence, dignity and choice in their lives.

Product development is driven by decades of experience with thousands of individual assessments bringing benefits to thousands of people. Contact us if you need a Neater Solution.

For our main product selection please see our website.

We offer **FREE HOME ASSESSMENTS** (UK) of our main equipment. This includes ensuring therapists and carers know how to make adjustments for optimum set up.

For our latest product information please contact us:



+44 1298 930123



info@neater.co.uk



www.neater.co.uk

| | |
|------------------------|----|
| Services | 04 |
| Neater Arm Support | 06 |
| NAS Additional | 10 |
| NAS Brackets | 13 |
| NAS-Zero | 15 |
| Neater Clamps | 17 |
| Neater Eater Manual | 18 |
| Neater Eater Robotic | 23 |
| Neater Crockery | 28 |
| Neater Cutlery | 30 |
| Neater Powered Drinker | 32 |
| Neater Drinkers | 34 |
| Neater Switches | 38 |
| Neater Holder | 40 |
| Neater Table | 45 |

FREE SERVICES

FREE HOME ASSESSMENTS (UK) of our main equipment. This includes ensuring therapists and carers know how to make adjustments for optimum set up. Occasionally we can facilitate short term trials.

3 Year Warranty included as standard. Return to base warranty - from date of invoice. UK Only. Excludes normal consumables e.g. switches, cutlery, crockery.

Technical Support as well as complete instruction manuals and instructional videos we offer technical support (UK) via your demonstrator or our main office, who will guide fixes via email, phone call or video call and if necessary advise on and arrange repairs.

FREE THERAPIST TRAINING

PROGRAMMES provide an overview of equipment and its suitability.

An experienced trainer will come to your chosen site and provide the following learning outcomes:

- Gain a clear overview of Neater equipment and features
- Learn to identify who would benefit from the equipment
- Technical Certificate
- Short sessions available to suit lunchtime meetings



College of
Occupational
Therapists
Approved

CHARGEABLE SERVICES

Visit: If you would like a Neater representative to visit and setup your Neater Equipment it can be arranged for this fixed fee within the UK.

Rental We offer rental programmes on specific products.

5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg tablet, batteries, switches, cutlery, crockery.

Refurbishment We offer refurbishment programmes on specific products up to 7 years.


Repairs We will repair Neater products up to 10 years.

NEATER ARM SUPPORT

NAS-1R/L Neater Arm Support - single unit complete - includes one setup visit (UK only) Available for both right and left arms. Right pictured here.



NAS-2 Neater Arm Support - double bracket unit complete - includes one setup visit (UK only)



NAS-A Standard attaching kit of parts for simple fitting of NAS to tubular frame wheelchair. Part dimensions vary dependent on wheelchair and require trained fitter



NAS-SL Neater Arm Support Sling Large, right or left handed. Can be used with any Neater Arm Support either powered or Zero power (sling frame, strap and fleece cover not included)



NAS-SS Neater Arm Support Sling Small, right or left handed. Can be used with any Neater Arm Support either powered or Zero power (sling frame, strap and fleece cover not included)



NAS-F Neater Arm Support Sling Frame, large or small. Can be used with any Neater Arm Support either powered or Zero power. Specify length required (sling, strap and fleece cover not included)



NAS-F Neater Arm Support Sling Small, right or left handed. Can be used with any Neater Arm Support either powered or Zero power (sling frame, strap and fleece



NAS-SF Sling Fleece Cover, right or left handed, large or small



NAS-VS Velcro strap (single)



NAS-TC Tube connector (pair)



...CONTINUED...

...CONTINUED...

NAS-Cable Power cable for Neater Arm Support



NAS-Pin — pin to join linked arm to NAS column



NAS-Arm Neater Arm Support linked arm complete. Available for both right and left arms



NASA-R/L-IS Neater Arm Support Arm, curved IN Top Arm, S shaped Mid Arm. Available for both right and left arms.



...CONTINUED...

NASA-R/L-IF Neater Arm Support Arm, curved IN Top Arm, Flat Mid Arm. Available for both right and left arms.



NASA-R/L-OS Neater Arm Support Arm, curved OUT Top Arm, S shaped Mid Arm. Available for both right and left arms.



NASA-R/L-OF Neater Arm Support Arm, curved OUT Top Arm, Flat Mid Arm. Available for both right and left arms.



NAS ADDITIONAL OPTIONS

NAS-Q Quick release system for NAS—allows quick removal/refit/transfer of complete NAS from wheelchair or frame/stand

NAS-Q2 Quick release system for NAS—including 2 holders

NAS-Assist Disengageable tilting mechanism for NAS to assist bringing hand to user's mouth automatically as the arm is lifted

NAS-P Separate stand for NAS

NAS-PC Separate stand for NAS with Casters



...CONTINUED...

NAS-Battery Neater Arm Support 24V battery and charger complete



NAS-EP 24V, 5A power supply for NAS



NAS-VS old Velcro strap (single) Standard size : 500mm



NAS-AA Adjustable angle fitting block—allows yoke bracket to be vertical when column is tilted. Available to fit right or left hand side



NAS-ET Extra tube for fitting – *Specify length*



NEATER ARM SUPPORT SERVICES

NAS-ReFit Re-fitting NAS to new Wheelchair or stand—UK only (assumes brackets reused)

NAS-5Year NAS 5 Year Package—includes: extension of warranty to 5 years and two additional fittings to new wheelchairs with any standard NAS brackets required. Conditions apply, UK only.

NAS-5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg batteries, switches.

NAS-Refurb Refurbish NAS: initial collection, full refurb to modern specification, electronics etc. (dependant on initial condition), new linked arms, 1 refit to wheelchair. 2 year warranty. (Cannot use 5 year plan).



NAS 5 Year Package

Overview

The **Neater Arm Support 5 Year Package (NAS-5Year)** has been developed in conjunction with clinicians in order to accommodate the changing needs of people using the Neater Arm Support - particularly young people.

The Neater Arm Support 5-year package includes:

- Extended warranty to 5 years.
- Up to two additional fittings of the Neater Arm Support within 5 years of the first fitting.
- Any additional standard brackets required to fit up to two new tubular frame wheelchairs.

The principal benefit of the 5-year package is that no additional funding needs to be sought within the period to enable fitting of the Neater Arm Support to new wheelchairs. Please contact us before purchasing new wheelchairs to ensure compatibility with the NAS.

One major advantage is that the user can have the Neater Arm Support fitted to their current chair even if they are expecting delivery of a new chair as it is important for the user to have the support early to help maintain movement in their arms. This can be compensated if the supply of a new wheelchair is very delayed. (Please note: a few weeks' notice may be required to make a date for fitting the Neater Arm Support.)

Terms and Conditions

All other Neater Arm Support terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. switches. Warranty excludes malfunctions from tampering and incorrect usage. Refit assumes re-use of original NAS fitting kit. Specialist bracketry required is chargeable. Refits are for original user only.



NAS 5 Year Warranty

Overview

The **Neater Arm Support 5 Year Warranty (NAS-5YW)** extends the warranty period to 5 years for your piece of mind.

Terms and Conditions

All other Neater Arm Support terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. switches. Warranty excludes malfunctions from tampering and incorrect usage. Warranty is return to base.



NAS Refurbishment

Overview

The **Neater Arm Support Refurbishment (NAS-Refurb)** is a full refurbishment to modern specification.

Includes:

- Initial collection.
- Refurbishment of electronics.
- New linked arms (dependant on initial condition).
- 1 refit to wheelchair.
- New 2 year warranty from date of delivery.

Terms and Conditions

All other Neater Arm Support terms and conditions apply. Read the safety instructions. UK Only. Warranty excludes consumables e.g. switches. Refit assumes re-use of original NAS fitting kit. Specialist bracketry required is chargeable. Refits are for original user only.



NAS FITTING KIT BRACKETS

NAS-CB-M5-22 Tube clamp body with M5 holes (standard) with 22mm straps




NAS-CB-M5-25 Tube clamp body with M5 holes (standard) with 25mm straps



NAS-CB-SLT Slotted tube clamp body – Comes with 4 tube connectors plus 2 additional 300mm tubes– Complex fitting



NAS-CB-P1 Permobil I type tube clamp body with wedges. Rarely used.



NAS-CB-6 Standard tube clamp body with counter bored 6mm through hole






NAS-CB-P40 CB-6 tube clamp body with small brass section for Puma 40 type chair

NAS-CB-P2 CB-6 tube clamp body with large block and brass sections for Permobil II type chair inc parts NAS-COR-EVOL-001 and NAS-COR-EVOL-002002

NAS-CB-M6-22 Tube clamp body with M6 end holes with 22mm straps

NAS-CB-M6-25 Tube clamp body with M6 end holes with 25mm straps

...CONTINUED



NEATER ARM SUPPORT ZERO POWER

NAS-Zero-R/L Zero-Power Neater Arm Support Right/Left Handed complete including table clamp kit. Specify sling & frame size

NAS0-HK Chair mounted parking hook for NAS0. Includes: NH-K, NH-008, 2X NH-E, NH-H, NH-C AND 1/4" UNC spacer.

NAS0-A1 Single arm NAS-0 zero powered arm support attaching kit. Includes 3 tube connectors, 1 NAS-252 docking connector, NAS-CB-M5-25, 2 attaching tubes, 2 plain tube 300 long.

NAS0-A2 Bilateral arm NAS-0 zero powered arm support attaching kit. Includes 6 tube connectors, 2 NAS-252 docking connectors, NAS-CB-M5-25, 2 attaching tubes, 4 plain tubes 300 long, crossbar.



NAS0-TC Neater Arm Support 0 Table Clamp Kit. Right or left handed to be specified.

NAS0-TCB Tube connecting base for NAS0. Can be used to mount the NAS0 on a separate stand tube.

NAS0-RENTAL Neater Arm Support 0 rental up to 1 year max. Customer to complete rental agreement. Full NAS0 value due in addition to rental charge if not returned. T&C's apply. UK only.

NAS0-5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg sling.

...CONTINUED...



Neater Arm Support Zero Rental

Overview

The **Neater Arm Support Zero Power Rental Package (NAS0-RENTAL)** has been developed in conjunction with clinicians in order to accommodate the changing needs of people using the Neater Arm Support Zero Power. Rental period can be arranged to a maximum of 1 year. Right or Left handed available.



Terms and Conditions
All other NAS0 terms and conditions apply. Read the safety instructions. UK Only. Warranty excludes consumables e.g. slings. Warranty excludes malfunctions from tampering and incorrect usage. Customer must complete rental agreement. Full NAS0 value due in addition to rental charge if not returned.



Neater.co.uk
Independence | Dignity | Choice

NAS0 5 Year Warranty

Overview

The **Neater Arm Support Zero Power 5 Year Warranty (NAS0-5YW)** extends the warranty period to 5 years for your piece of mind.

Terms and Conditions

All other NAS0 terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. slings. Warranty excludes malfunctions from tampering and incorrect usage. Warranty is return to base.



Neater.co.uk
Independence | Dignity | Choice

CLAMPS

NES-C Clamp (without Large Foot) (single) these come with NEM as standard. Available in black or white.




NES-CS Small Clamp (without Large Foot) (single). White only.



NES-CD Large Clamp for deeper lip tables (up to 110mm) (without Large Foot) (single). White only.



NES-CTB Replacement "tommy" bar for the Neater Clamp. Includes ball knob ends.



NEM MODULAR OPTIONS

NEM-R/L Multi adjustable damped arm with forward & back stops, slimline board, spoon lift Neater Eater, Right and left handed options available

NEM-C-R/L Neater Eater Manual , includes cutlery & crockery. Right and left handed options available. Right is shown here,

NEM-T-R/L Neater Eater Manual, Shorter Arm. Right and left handed options available

NEM-5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg crockery & cutlery.

NEM-Refurb Full workshop service, new cutlery & crockery and return via courier. New belt and spring, and replace any parts to meet modern standards. Extends the warranty by 12 months.



NEM 5 Year Warranty

Overview

The **Neater Eater Manual 5 Year Warranty (NEM-5YW)** extends the warranty period to 5 years for your piece of mind.

Terms and Conditions

All other NEU terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. crockery and cutlery. Warranty excludes malfunctions from tampering and incorrect usage. Warranty is return to base.



NEM Refurbishment

Overview

The **Neater Eater Manual Refurbishment (NEM-REFURB)** is a full refurbishment to modern specification.

Includes:

- Initial collection.
- Full workshop service.
- New belt and spring.
- New cutlery & crockery.
- New 2 year warranty from date of delivery.



Terms and Conditions

All other Neater Eater Manual terms and conditions apply. Read the safety instructions. UK Only. Warranty excludes consumables e.g. crockery & cutlery. Warranty excludes malfunctions from tampering and incorrect usage.



Parts below can be ordered within NEM code (eg NEM-GH2JPS-L) or separately with prefix NEMA-

NEMA-C Standard cutlery and crockery package. Includes 1 of each: NEP-L, NEP-S, NEP-R, NEC-FL, NEC-FS.

NEMA-C2/3 Cutlery package including any 2 /3 plates, any 2 /3 cutlery and a plate stand. UK only, excludes partitioned plate.

NEMA-S-R/L Forward Spring attachment for Right/Left handed

NEMA-J 40mm Baseboard raising spaces (reduces plate to mouth distance) single unit

...CONTINUED...



...CONTINUED...

NEMA-F Forward mounting baseboard extension (to bring arm or plate forward)



NEMA-G Plate flick stop guard



NEMA-H Adjustable Handle (quick release)



NEMA-P Remote plate Turner



NEMA-CP Colour Pink



Parts below are ordered as spares with prefix NEMS-

NEMS-B Ultra slim narrow bare board flush fit



NEMS-E Clamp offset board (to move clamp away from awkward positions) board replaces standard board



NEMS-BP Circular Baseboard Disk to fit under plate



NEMS-H Cover H



...CONTINUED...

NEMS-S-R/L Forward spring alone—
Right or left handed

NEMS-U Counter balance spring
alone

...CONTINUED



NEATER EATER ROBOTIC

NER Neater Eater Robotic complete. Chargeable crockery and cutlery.

NER-C Neater Eater Robotic - Electric Programmable Version complete incl. Plates, plate stand & cutlery. Comes with silver case as standard. Includes carry bag. Includes cutlery and crockery: NEP-L, NEP-S, NEP-R, 1x NEC-FL, 1x NEC-FS.

NER-P Neater Eater Robotic - Electric Programmable Version complete incl. Plates, plate stand & cutlery. Comes with pink case and carry bag. Chargeable cutlery and crockery.

NER-Rental Neater Eater Robotic rental any period to 1 year max. Customer to complete form to agree to appropriate return of equipment. If not returned, the full value of a new NER to be paid in addition to the rental charge.



NER-Rental



Overview

The Neater Eater Robotic Rental Package (NER-Rental) has been developed in conjunction with clinicians in order to accommodate the changing needs of people using the Neater Eater Robotic.

Terms and Conditions

All other Neater Eater Robotic terms and conditions apply. Read the safety instructions. Customer to complete form to agree to appropriate return of equipment. If not returned the full value of a new NER to be paid in addition to the rental charge. Rental period can be arranged to a maximum of 1 year.



Neater.co.uk
Independence • Dignity • Choice

NER-R2P Purchasing of a rental unit after the rental period expires. Technical support or return to Neater as required in order to ensure the device is up to modern standard. Includes 3 year warranty from delivery date.

NER-5MP 5 Year Maintenance Plan. Includes extended return to base 5 year warranty (from date of invoice) and service between years 3 and 4. Must be purchased with new unit. UK Only. Excludes consumables eg tablet, switches, crockery & cutlery. Must be purchased with new unit. UK Only.

NER-Refurb Refurbishment of Neater Eater Robotic up to meet modern standard. Includes new cutlery, crockery, batteries, tablet and cover. Belt tightening, calibration, software updates. Comes with 2 year warranty.

NER-C2/3 Neater Eater Robotic - Electric Programmable Version complete incl. Plates, plate stand & cutlery. Comes with silver case as standard. Includes carry bag. Includes ANY 2/3x cutlery and crockery plus 1x plate stand. Excludes partitioned plate. UK ONLY.

...CONTINUED...

NER-Rent to Purchase

Overview

Purchasing of a **Neater Eater Robotic** rental unit after the rental period expires. Technical support or return to Neater as required in order to ensure the device is up to modern standard. Includes 3 year warranty from original delivery date.

Terms and Conditions
All other Neater Eater Robotic terms and conditions apply. Read the safety instructions. Warranty excludes consumables e.g. tablet, switches, crockery & cutlery. Warranty excludes malfunctions from tampering and incorrect usage.



5 Year Maintenance Plan

Overview

The **Neater Eater Robotic** 5 Year Maintenance Plan. Includes extended return to base 5 year warranty (from date of invoice) and service.



Contact us to arrange your service. We will send a loan unit to you while your robot is returned to us. We will check everything inside the robot is working correctly, recalibrate and replace parts as necessary. Your robot will then be returned to you, with the confidence it will continue to work well for years to come.

Terms and Conditions
All other Neater Eater Robotic terms and conditions apply. Read the safety instructions. Your service needs to be arranged for between 9am and 5pm, 7 days a week. Can only be purchased with new unit. UK Only. Warranty and service excludes consumables e.g. tablet, switches, crockery & cutlery. Warranty and service excludes malfunctions from tampering and incorrect usage.



NER-Refurb

Overview

Refurbishment of **Neater Eater Robotic** up to meet modern standard. Includes new cutlery, crockery, batteries, tablet, cover, belt tightening, calibration and software updates. Comes with 2 year warranty.

Terms and Conditions
All other Neater Eater Robotic terms and conditions apply. Read the safety instructions. Warranty excludes consumables e.g. tablet, switches, crockery & cutlery. Warranty excludes malfunctions from tampering and incorrect usage.



NER-V Replacement covers for the Neater Eater Robotic - Standard Metallic Silver.



NER-VP Replacement covers for the Neater Eater Robotic -Metallic Pink.



NER-VT Replacement Top Arm Cover for the Neater Eater Robotic



NER-VM Replacement Mid Arm Cover for the Neater Eater Robotic



NER-VL Replacement Low Arm Cover for the Neater Eater Robotic



...CONTINUED...

...CONTINUED...

NER-VB Replacement Main Body Cover for the Neater Eater Robotic



NER-EP 12v regulated power supply for Neater Eater Robotic. Includes low voltage transformer and mains power lead.



NER-BAG Bag for NER with custom foam inserts for additional product protection



NER-CASE Black plastic protective case for transporting NER, comes with custom foam cutout to fit around the robot for protection



NER-Tablet Replacement tablet for Neater Eater Robotic



...CONTINUED

NER-BATTERY Replacement battery pack for the Neater Eater Robotic. Requires professional installation.



NER-PH Replacement plate holder for Neater Eater Robotic



NEATER PLATES

All plates have spigots

NEP-S Curled-over deep lip ceramic plate with spigot - Small 21 x 3.3cm deep

NEP-L Curled-over deep lip ceramic plate with spigot - Large 26 x 2.5cm deep

NEP-P Plastic plate with vertical lip with spigot. Blue or White (we do not recommend putting plastic plates in the dishwasher as this will shorten their lifespan.)

NEP-B Soup Bowl with saucer attached with spigot



...CONTINUED...

NEP-R Plate stand (ceramic)



NEP-PR Plate stand (plastic)



NEP-H Heated Ceramic Plate



NEP-P Partitioned plate . To work with the NER you must order programmable code.



NEATER CUTLERY

NEC-FL Large flat front co-polymer silver Neater Spoon (standard)



NEC-FS Small flat front co-polymer silver Neater Spoon (standard)



NEC-ML Flat front stainless steel spoon



NEC-MSN Narrow flat front stainless steel spoon



NEC-MSP Soup Spoon made from stainless steel



...CONTINUED...

NEC-MK "Spork" (spoon with fork prongs) made with stainless steel



NEC-C Small Round Fronted Maroon Spoon



NEC-MXL Extra Large Flat Front Stainless Steel Spoon



NEC-H Cutlery Holder alone



NEB-S2 Sandwich & Snack Holder



NEATER POWERED DRINKERS

NPD-D Neater Powered Drinker including peristaltic pump, 2 flexible straws, straw holder, switch and battery charger

NPD-5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg batteries, straw, bag, switch.

NSM Electric Straw Mover includes bag, 2 x flexible straws, non-return valve, straw cleaning brush, charger.

NSM-5YW 5 Year extended return to base warranty - extends the warranty period to 5 years (from date of invoice). Must be purchased with new unit. UK Only. Excludes normal consumables eg batteries, switch.

NSM-A Electric Straw Mover with Neater Holder Arm (instead of bendy arm) includes bag, 2 x flexible straws, non-return valve, straw cleaning brush, charger (UK only).



NPD 5 Year Warranty

Overview

The Neater Powered Drinker 5 Year Warranty (NPD-5YW) extends the warranty period to 5 years for your piece of mind.

Terms and Conditions

All other NPD terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. straws and switches. Warranty excludes malfunctions from tampering and incorrect usage. Warranty is return to base.



NSM 5 Year Warranty

Overview

The Neater Straw Mover 5 Year Warranty (NSM-5YW) extends the warranty period to 5 years for your piece of mind.

Terms and Conditions

All other NSM terms and conditions apply. Read the safety instructions. Can only be purchased with new unit. UK Only. Warranty excludes consumables e.g. straws and switches. Warranty excludes malfunctions from tampering and incorrect usage. Warranty is return to base.



NEATER POWERED DRINKER SPARES



NPD-S/4S 4x Spare flexible straws for Neater Powered Drinker. Each straw 120cm long (precision platinum cured).

ND-SB Long Neater Straw (120cm long with non-return valve)

NPD-25S 25 meter roll of flexible straw material used with the Neater Powered Drinker, Neater Straw Mover and Neater Steady Straw (plus other products as required).

NPD-B Spare rechargeable battery pack for Neater Powered Drinker.

NPD-BAG Bag for Neater Powered Drinker 27x19x14cm or straw mover.



NEATER DRINKER OPTIONS & SPARES

ND-D Neater Drinker including ceramic mug, straw, non-return valve, mouthpiece & syringe. Available in blue, white, green or black



ND-T Neater Drinker Thermo Capacity: 600ml. (1 pint)



ND-M Neater Drinker Metallic Capacity: 400ml. Fits ND-MB stable base and ND-MH cup holder—attaches to wheelchairs etc.



ND-KV2 Neater Steady-Double-Straw for use with any cup — including 2 flexible straws, bendable support with table clamp, 2 non-return valve and FREE straw cleaning brush



ND-KVC2 Neater Steady-Double-Straw Tube Clamp version- for use with any cup - including 2 flexible straws, bendable support (NOT Neater Holder Arm) with c-clamp, ball head, 2 non-return valves, 2 spare straw clips and FREE Straw Cleaning Brush.



...CONTINUED...

ND-MKV Neater Cup Holder with Steady straw




ND-MKVC Neater Cup Holder with Steady-Straw - Tube Clamp version.



ND-MH Neater Cup Holder



ND-SC Straw Cleaning Brush



ND-KS Spare flexible for Neater Steady-straw 80cm long



...CONTINUED...

ND-S Neater Straw (60cm long with non-return valve)




ND-SB Long Neater Straw (120cm long with non-return valve)




ND-SW 60cm Straw alone without non-return valve



ND-SBW 120cm Straw alone without non-return valve



ND-SV Spare Non-Return Valve for Straw



...CONTINUED

ND-C Neater Straw Clip



ND-CS 120mm (12cm) length of hard straw without non return valve



NEATER SWITCHES

NB-N Neater Switch




NB-F Foot Switch B218



NB-N3H PUSH BUTTON (Neater Switch) plate for holding up to 3 Neater Switches - includes non-slip feet (use Velcro included with Neater Switches to attach) (12cm x 5cm) (Neater Switches ONLY, not for use with Foot Switches)




NB-F3H FOOT SWITCH plate for holding up to 3 Foot Switches - includes non-slip feet and strong Velcro (24cm x 10cm) (Foot Switches ONLY, not for use with Neater Switches)



NB-N3 Switch pack of 3x Neater Switches to use as required. Optional choice of colour subject to availability.



...CONTINUED



NB-V Spare sticky-back base two-part Velcro ring

NEATER HOLDER

The Neater Holder is a modular system for firmly holding objects at adjustable positions and angles.

It must be ordered as a complete system.

The Universal Holder is suitable for objects such as an electric razor or special all-angle electric toothbrush.

To order your Neater Holder it's as simple as 1..2..3!

1. Choose the base
2. Choose the body
3. Choose the holder (s).

1. Choose your base

NH-B Neater Holder Base - requires Neater Clamp (small, standard or large, not included)

NH-C Neater Holder C-Clamp - for clamping to tubes 1cm for 4cm diameter

NH-W Wall mount - adjustable angle with screws to fix to wall (by suitably skilled person)



2. Choose the Body

NH-A Neater Holder Arm - angle adjustable sturdy arm requires base, clamp, extensions and attachments (not included) 2x 11cm arms with ball adjust each end and locking knob.

Alternatively for a more compact, but less adjustable with less range:

NH-H Ball head - adjustable angle (probably not needed if using Neater Holder Arm)

You may also want with either of the above, or straight into any base:

NH-E Neater Holder Extension - 10cm long extension piece

3. Choose your Holder (s)

NH-U Neater Holder Universal Clamp - clamps objects up to 9cm diameter (*CHECK*)


NH-T Neater Tool Holder Universal Clamp - fits snack holder (not included) and other Neater cutlery

...CONTINUED...



...CONTINUED...

NH-S1 Single Switch Holder



NH-S2 2/3 Switch Holder



NH-P1 Phone holder 5.8cm to 10cm width - to fit Neater Holder system



NH-P2 Tablet holder



NH-NERT NER Tablet Holder - with magnets that hold Neater Eater Robotic tablet



NH-M Bracket for cupholder and straw mover to fit Neater Holder system - Fix cupholder and/or Neater Straw Mover to Neater Holder Arm or C-Clamp



NH-V2 Double Straw Holder and Accessories for fitting to Neater Holder system - for use with any cup - including 2 flexible straws, 2 non-return valve and FREE Straw Cleaning Brush. Requires clamps and Neater Holder arm etc (not included)

NH-K Parking Hook for NAS-0 - for use when chair mounting (rather than table mounting) NAS-0.

NH-NPDH Neater Powered Drinker Holder attachment for the Neater Holder System. Includes NH-007A NPD Holder, NH-009 NPD, NH-C x2, NH-H x2, NH-E x2. To be used with the Neater Holder System.

...CONTINUED...



...CONTINUED...

Examples

It's helpful when ordering to list the parts starting at the base so that the product can be configured correctly for delivery.

Table mounted snack holder:

1. BASE: NC-C (Black); NH-B;
2. BODY: NH-A; NH-E;
3. HOLDER: NH-T; NEB-S2



Tube mounted universal holder

1. BASE: NH-C;
2. BODY: NH-A; NH-E;
3. HOLDER: NH-U;




Table Mounted NER Magnetic Holder

1. BASE: NC-C (Black); NH-B;
2. BODY: NH-A;
3. HOLDER: NH-NERT



Cup Holder & Straw Mover with chair mounting

NSM & ND-MH Straw Mover & Cup Holder

1. BASE: NH-C;
2. BODY: NH-H
3. HOLDER: NH-M



Chair mounted Parking Hook (for NAS-0)

Usually this would be on a wheelchair. Can be ordered as **NASO-HK**

1. BASE: NH-C;
2. BODY: NH-H, NH-E (x2)
3. HOLDER: NH-K



NEATER TABLE

NT-TC Standard Neater Table



...CONTINUED



Neater Solutions Ltd

12 Burlington Road,

Buxton,

Derbyshire,

SK17 9AL

UK

Tel: **+ 44 1298 238823**

Email: **info@neater.co.uk**

www.neater.co.uk

221018

