

Order Form

PROGEO UK (English)

Company stamp

TEKNA ADVANCE

Version WITH FOOTRESTS REMOVABLE AND SWING AWAY

Date

User name

Responsible

Folding frame wheelchair

Wheelchair measures

Seat width (LS) ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48

Seat depth (PS) ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. - calculated with angle backrest 90° to the seat

Footplate distance (DP) (Real distance between seat and footplate)

IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)

Setting ☒ Standard ☐ Prudential ☐ Active

Camber ☒ 2° ☐ 0° ☐ 4°


Backrest Height (TS) Std.Value 33,5(min 24.5-max 47)

Backrest angle (to the seat) ☒ 90° ☐ -6° ☐ +6°


Front height (STD 49) - Rear height (STD 43)

Front frame


Frame Angle



☒ 110°




☐ 120°



☐ 95°

Rear frame

☒ Aluminium standard frame



☐ Short frame - Ideal with Seat depth (PS)37,5 - 40

☐ Long frame - Ideal with Seat depth (PS)42,5 - 45 - 47,5

Frame colour

☐ Bright silver

☐ Opaque black

☐ Pearl white

Backrest

☒ Handle tube STD aluminium

☐ Backrest with folding handles



02958

£ 218

Casters

☒ 4" Soft

☐ 3" Sport

☐ 5" Soft

☐ 6" Soft

☐ 6" Pneumatic

Not available with frame angle 95°

Not available with frame angle 95°

☐ 3" Active Soft (wide)

(Aluminium alloy rim)

06824

£ 109

☐ 4" Active Sport (narrow)

(Aluminium alloy rim)

01398

£ 109

☐ 4" Active Soft (wide)

(Aluminium alloy rim)

03943

£ 109

☐ 5" Active Soft (wide)

(Aluminium alloy rim)

07518

£ 109

Support front fork

☒ Standard support front fork



Colour set: Rear wheels and forks - To be matched only with:Rear wheels Ultra-light Progeo - Standard! (22 - 24)

☒ Polished black 22" 24" 25" 26"

☐ Bright silver 22" 24"

00833

£ 109

☐ Bright anthracite 24"

00833

£ 109

☐ Bright orange 22" 24"

00833

£ 109

☐ Bright fuchsia 22" 24"

00833

£ 109

☐ Bright blue 22" 24"

00833

£ 109

☐ Bright red 22" 24"

00833

£ 109

☐ Bright lime 22" 24"

00833

£ 109

☐ Pearl white 24"

00833

£ 109

English
Date Order User name Address

Rear wheels

- ☆ 24" Ultra-light Progeo - Standard!
- ☐ 22" Ultra-light Progeo - Standard!
- ☐ 25" Ultra-light Progeo - Standard!
- ☐ 26" Progeo - Standard!
- ☐ 24" ULTRA-LIGHT PROGEO KEY SPORT **Weight reduction approx. 450gr.** 04383 **£ 742**
- ☐ 24" Ultra-light "Spinergy LX" 01685 **£ 742**
- ☐ 24" Ultra-light "Spinergy Flex rim" 02270 **£ 1200**
- ☐ 24" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes ☐ Hub and blade spokes black 07846 **£ 742**
- ☐ 25" Ultra-light "Spinergy LX" 02101 **£ 742**
- ☐ 25" Ultra-light "Spinergy LXL" ☐ Black hub and white blade spokes 07847 **£ 742**
- ☐ 26" Ultra-light "Spinergy LX" 02102 **£ 742**
- ☐ 24" Off-road (aluminium handrim black) 01392 **£ 164**
- ☐ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°
- ☐ Right ☐ Left 02103 **£ 327**

Rear wheel tyre

- ☆ High pressure tyre (22"-24"-25"-26")
- ☐ High pressure profile tyre (24")
- ☐ Marathon Plus (24"-25") **with anti-puncture protection** 01669 **£ 55**
- ☐ Solid tyre Kik Mako - Colour Black (24") 01041 **£ 87**
- ☐ Solid tyre Kik Mako - Colour Grey (24") 01751 **£ 87**

Handrim - Wheel distance

☆ 3 cm ☐ 4 cm

- ☆ Aluminium (22"-24"-25"-26") 09170 **£ 109**
- ☐ Aluminium black SOFT-TOUCH (24"-25") 01684 **£ 109**
- ☐ Chromed Steel (24"-25") 00016 **£ 218**
- ☐ Titanium (24"-25") 00017 **£ 262**
- ☐ Ultra Grepp (24"-25") 00834 **£ 305**
- ☐ Max Grepp (24"-25") 00632 **£ 305**
- ☐ Ergo Para (24") 06427 **£ 316**
- ☐ SURGE LT (24"-25") **Oval countour with a rubber insert** 00693 **£ 175**
- ☐ Quad (24") 04050 **£ 142**
- ☐ Silicon cover for standard aluminium handrim Available for measures (22"-24"-25"-26")

Wheel locks

- ☆ Bent lever push to lock
- ☐ Straight lever Push to lock
- ☐ Straight lever Pull to lock
- ☐ Sport brake OUT-FRONT
- ☐ Brake with extension lever push to lock 01949 **£ 55**
- ☐ Brake with extension lever pull to lock 03372 **£ 55**

Footplate - See reference tables N°3


- ☆ Plastic footrest rake adjustable
- Minimum footplate distance 35cm.** ☆2/3 Internal ☐ 2/3 External
- ☐ One piece footplate - Flip-up! rake and depth adjustable
- Minimum footplate distance 35cm.** ☐ Internal ☆2/3 Internal ☐ External
- ☐ **Minimum footplate distance 30cm. Available only with NEW!!!Multi-function and rotative system - 03945 £ 164,00**
- ☐ One piece aluminium footplate - Flip-up! rake and depth adjustable
- Minimum footplate distance 35cm.** ☐ Internal ☆2/3 Internal ☐ 2/3 External ☐ External
- ☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up.
- Minimum footplate distance 35cm.** ☐ Internal ☐ External 03291 **£ 240**

Date Order English User name Address

Mudguard

Mudguard

- ☒ ABS mudguard
- ☐ ABS Straight sideguard
- ☐ Carbon fibre mudguard
- ☐ Carbon fibre Straight sideguard



03946 £ 262
08918 £ 262

Removable version

- ☐ Removable ABS mudguard
- ☐ Removable ABS fibre straight sideguard
- ☐ Removable carbon fibre mudguard
- ☐ Removable Carbon fibre straight sideguard



08921 £ 109
08921 £ 109
08919 £ 371
08920 £ 371

Armrests

- ☐ Swing away armrests
- ☐ Removable armrest with ABS straight sideguard
- ☐ Removable armrest with carbon-fiber straight sideguard



00694 £ 218



08881 £ 196
08882 £ 393

Kind of pad to be mounted on chosen armrest

- ☒ Pad STANDARD
- ☐ Ergonomic pad (30x6 cm.)

04104 £ 55

Optional 1

Only one optional devices for side can be chosen

- ☐ Transit wheels Exclude the installation of :
Anti-tippers -Tipping aid - Crutches support
- ☐ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance)
☐ Right
☐ Left
- ☐ Tipping aid ☐ Right
☐ Left
- ☐ Crutches support ☐ Right
☐ Left



00836 £ 131



00023 £ 109
00024 £ 109



00025 £ 55
01046 £ 55



00838 £ 65
00837 £ 65

Optional 2

- ☐ Tetraclip (Quick release axle for quadriplegic)
- ☐ Quick release front forks To be matched only with: Standard support front fork
- ☐ Pad for front frame protectors
- ☐ Push handles for companion (Removable and height adjustable)
- ☐ Safety pelvic belt
- ☐ Harness
- ☐ Lumbar support (to insert in the backrest upholstery)
- ☐ Swing away lateral support ☐ Right
☐ Left
Not available with Posture backrest PHYSIO
- ☐ Spoke guards Black (22"-24")
- ☐ Spoke guards Transparent (22"-24")
- ☐ Spoke guards Progeo design (22"-24")
- ☐ Spoke guards Dragon (22"-24")



01048 £ 153

03972 £ 55

04103 £ 109

00013 £ 33



02097 £ 218
02098 £ 218

01043 £ 131

00835 £ 131

00019 £ 153

02106 £ 153

Additional Kit

- ☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour)
- ☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour)

03378 £ 436

09357 £ 764

Date Order English User name Address

PHYSIO - Posture backrest

Std. Backrest upholstery not provided

☐ Aluminum shell - Lateral retaining STANDARD 7,5cm Height ☐ 25 ☐ 33 ☐ 41 05308 £ 884



Height ☐ 52 05309 £ 938

Lateral retaining wing ☐ A=high ☐ B=low

☐ Aluminum shell - Lateral retaining DEEP 15cm Height ☐ 25 ☐ 33 ☐ 41 05310 £ 938



H25 not available for seat width 39-42-45-48 05311 £ 993

Height ☐ 52

Lateral retaining wing ☐ A=high ☐ B=low

☐ Carbon shell - Lateral retaining STANDARD 7,5cm Height ☐ 25 ☐ 33 ☐ 41 05312 £ 1015



Lateral retaining wing ☐ A=high ☐ B=low


COVER

☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




☐ SWING AWAY - Lateral support
Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width,
adjustable (0 ÷ -1,5cm)




Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647
Small (14x9,5)
☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width,
adjustable (-4 ÷ -5,5cm)



Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647
Small (14x9,5)
☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width,
adjustable (-8,5 ÷ -10cm)



Large (17,5x11,5)
☐ Left 05646 ☐ Right 05647 £ 175
Small (14x9,5)
☐ Left 05642 ☐ Right 05643 £ 175
(with additional kit)

☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

☐ Headrest adjustable - Double ball mounting system

Not available for Physio backrest's height 25 and with carbon shell






☐ Large (34x14) 05648 £ 393
☐ Small (28x14)

Date Order English User name Address

Note

Enclosure to wheelchair order form Tekna Advance - Version with footrests removable and Swing away

REFERENCE TABLES

| 1 | | Front height | | | |
|---|--|-------------------|-------------------|-------------------|-------------------|
| With closed frame (95°) max 5" caster | | | | | |
| FORK | | Casters | | | |
| | | 3" | 4" | 5" | 6" |
| SMALL  | | min. 39 max 42 | min. 40 max 43 | min. 42 max 43 | min. 44 max 45 |
| MEDIUM  | | min. 43 max 45 | min. 44 max 46 | min. 44 max 47 | min. 46 max 48 |
| LONG  | | min. 46 max 52 | min. 47 max 53 | min. 46 max 54 | min. 47 max 55 |
| The size of the fork will be defined only by Rehateam Progeo in accordance with wheelchair final configuration. | | | | | |

| 2 | Rear height | | | |
|-------|-------------|---------|---------|---------|
| | Rear wheels | | | |
| | 22" | 24" | 25" | 26" |
| H min | min. 36 | min. 38 | min. 39 | min. 40 |
| H max | max 44 | max 48 | max 49 | max 51 |

| 3 | Footplate option | Footplate distance |
|--|--|--------------------|
| | Plastic footrest | min. 35 |
| | One piece footplate - Flip-up!, with multifunction rotative system | min. 30 |
| | One piece aluminium footplate - Flip-up! | min. 35 |
| | One piece footplate - Flip-up! | min. 35 |
| | One piece Carbon fibre footplate | min. 35 |
| Calculated with the best configuration and with the footplates positioned on inside. | | |
| IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.) | | |

| 4 | Rear wheels | |
|---|-------------------------------------|--|
| | Model | Features |
| | 22" Ultra-light Progeo - Standard! | (not available in anthracite color) |
| | 24" Ultra-light Progeo - Standard! | (available for all coulores) |
| | 24" ULTRA-LIGHT PROGEO KEY SPORT | (silver hub, black rim and 16 spokes) |
| | 24" Ultra-light "Spinergy LX" | (silver hub, black rim and spokes) |
| | 24" Ultra-light "Spinergy LXL" | (black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes) |
| | 24" Ultra-light "Spinergy Flex rim" | (silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre) |
| | 24" Off-road | (colour: black) - (aluminium handrim - tyre Off-road) |
| | 25" Ultra-light "Spinergy LX" | (silver hub, black rim and spokes) |
| | 25" Ultra-light "Spinergy LXL" | (black hub and rim, 18 white blade spokes) |
| | 26" Progeo - Standard! | (colour: black) |
| | 26" Ultra-light "Spinergy LX" | (silver hub, black rim and 18 spokes) |