

Date .....

User name .....

Company stamp

Folding frame wheelchair

Responsible .....

**Wheelchair measures**

Seat width (LS) ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48

Seat depth (PS) ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5  
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. - calculated with angle backrest 90° to the seat

Footplate distance (DP) ..... (Real distance between seat and footplate)  
**IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)**

Setting ☒ Standard ☐ Prudential ☐ Active

camber ☒ 2° ☐ 0° ☐ 4°


Backrest Height (TS) Std.Value 33,5 .....(min 24.5-max 47)

Backrest angle (to the seat) ☒ 90° ☐ -6° ☐ +6°

Front height ..... (STD 49) - Rear height ..... (STD 43)


**Front frame**

☆ Aluminium standard frame Angle 100°



- ☐ Short frame - Ideal with Seat depth (PS)37,5 - 40
- ☐ Long frame - Ideal with Seat depth (PS)42,5 - 45 - 47,5

☐ Aluminium frame Angle 90° - (max 5" caster)



- ☐ Short frame - Ideal with Seat depth (PS)37,5 - 40
- ☐ Long frame - Ideal with Seat depth (PS)42,5 - 45 - 47,5

**Frame colour**

- ☐ Bright silver
- ☐ Opaque black
- ☐ Pearl white
- ☒ Opaque orange  
08341 £ 55

**Backrest**

☆ Handle tube STD aluminium

☐ Backrest with folding handles



02958 £ 218

**Casters**

☆ 4" Soft

- ☐ 3" Sport
- ☐ 5" Soft
- ☐ 6" Soft
- ☐ 6" Pneumatic

(Available only with 100° angle)  
(Available only with 100° angle)

- ☐ 3" Active Soft (wide) (Aluminium alloy rim) 06824 £ 109
- ☐ 4" Active Sport (narrow) (Aluminium alloy rim) 01398 £ 109
- ☐ 4" Active Soft (wide) (Aluminium alloy rim) 03943 £ 109
- ☐ 5" Active Soft (wide) (Aluminium alloy rim) 07518 £ 109

**Support front fork**

☆ Standard support front fork



☐ Front fork support V Design



06850 £ 131

Date	Order	English User name	Address
------	-------	----------------------	---------

**Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)**

☆ Polished black	22" 24" 25" 26"		
○ Bright silver	22" 24"	00833	£ 109
○ Bright anthracite	24"	00833	£ 109
○ Bright orange	22" 24"	00833	£ 109
○ Bright fuchsia	22" 24"	00833	£ 109
○ Bright blue	22" 24"	00833	£ 109
○ Bright red	22" 24"	00833	£ 109
○ Bright lime	22" 24"	00833	£ 109
○ Pearl white	24"	00833	£ 109

**Rear wheels**

☆ 24" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!	(colour: black)		
<input type="checkbox"/> 26" Progeo - Standard!	(colour: black)		
○ 24" ULTRA-LIGHT PROGEO KEY SPORT	Weight reduction approx. 450gr.	04383	£ 742
○ 24" Ultra-light "Spinergy LX"		01685	£ 742
○ 24" Ultra-light "Spinergy Flex rim"		02270	£ 1200
○ 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes <input type="checkbox"/> Hub and blade spokes black	07846	£ 742
○ 25" Ultra-light "Spinergy LX"		02101	£ 742
○ 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white blade spokes	07847	£ 742
○ 26" Ultra-light "Spinergy LX"		02102	£ 742
○ 24" Off-road (aluminium handrim black)		01392	£ 164
○ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°	<input type="radio"/> Right <input type="radio"/> Left	02103	£ 327

**Rear wheel tyre**


☆ High pressure tyre	(22"-24"-25"-26")		
<input type="checkbox"/> High pressure profile tyre	(24")		
○ Marathon Plus	(24"-25")	with anti-puncture protection	01669 £ 55
○ Solid tyre Kik Mako - Colour Black	(24")		01041 £ 87
○ Solid tyre Kik Mako - Colour Grey	(24")		01751 £ 87

**Handrim - Wheel distance**

 ☆ 3 cm    ☐ 4 cm

☆ Aluminium	(22"-24"-25"-26")		
○ Aluminium black SOFT-TOUCH	(24"-25")	09170	£ 109
○ Chromed Steel	(24"-25")	01684	£ 109
○ Titanium	(24"-25")	00016	£ 218
○ Ultra Grepp	(24"-25")	00017	£ 262
○ Max Grepp	(24"-25")	00834	£ 305
○ Ergo Para	(24")	00632	£ 305
○ SURGE LT	(24"-25")	Oval countour with a rubber insert	06427 £ 316
○ Quad (24")	(24")	00693	£ 175
○ Silicon cover for standard aluminium handrim	Available for measures (22"-24"-25"-26")	04050	£ 142

**Wheel locks**

☆ Bent lever push to lock			
<input type="checkbox"/> Straight lever Push to lock			
<input type="checkbox"/> Straight lever Pull to lock		06428	£ 109
○ Sport brake OUT-FRONT		01949	£ 55
○ Brake with extension lever push to lock		03372	£ 55
○ Brake with extension lever pull to lock			

Date ..... Order ..... English User name ..... Address .....

## Footplate - See reference tables N°3

☆ Automatic Closure footplate rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal ☆ 2/3 Internal ☐ External




☐ Minimum footplate distance 30cm. Available only with NEW!!!Multi-function and rotative system - 03945 £ 164,00

☐ One piece footplate - Flip-up! rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal ☆ 2/3 Internal ☐ External



☐ Minimum footplate distance 30cm. Available only with NEW!!!Multi-function and rotative system - 03945 £ 164,00

☐ One piece aluminium footplate - Flip-up! rake and depth adjustable

Minimum footplate distance 35cm.


☐ Internal ☆ 2/3 Internal ☐ 2/3 External ☐ External



☐ Plastic footrest rake adjustable

Minimum footplate distance 35cm.

☆ 2/3 Internal ☐ 2/3 External



06861

£ 76

Need a stretch calf band.

☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.

☐ Internal ☐ External



03291

£ 240

## Mudguard


### Mudguard

☆ ABS mudguard

☐ ABS Straight sideguard

☐ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard



03946

£ 262

08918

£ 262

### Removable version

☐ Removable ABS mudguard

☐ Removable ABS fibre straight sideguard

☐ Removable carbon fibre mudguard

☐ Removable Carbon fibre straight sideguard



08921

£ 109

08921

£ 109

08919

£ 371

08920

£ 371

## Armrests

☐ Swing away armrests

(Swing away, removable and height adjustable armrest)



00694

£ 218

☐ Removable armrest with ABS straight sideguard



08881

£ 196

☐ Removable armrest with carbon-fiber straight sideguard



08882

£ 393

## Kind of pad to be mounted on chosen armrest

☆ Pad STANDARD

☐ Ergonomic pad (30x6 cm.)

04104

£ 55

## Optional 1

### Only one optional devices for side can be chosen

☐ Transit wheels

Exclude the installation of :

Anti-tippers -Tipping aid - Crutches support



00836

£ 131

☐ Anti-tipper with rapid hookup

(inside rotation to reduce the encumbrance)

☐ Right

☐ Left



00023

£ 109

00024

£ 109

☐ Tipping aid

☐ Right

☐ Left



00025

£ 55

01046

£ 55

☐ Crutches support

☐ Right

☐ Left



00838

£ 65

00837

£ 65

Date ..... Order ..... English User name ..... Address .....

## Optional 2

<input type="checkbox"/> Tetraclip	(Quick release axle for quadriplegic)	03947	£ 131
<input type="checkbox"/> Quick release front forks	To be matched only with: Standard support front fork	00030	£ 131
<input type="checkbox"/> Pad for front frame protectors		07058	£ 39
<input type="checkbox"/> Push handles for companion	(Removable and height adjustable)	01048	£ 153
<input type="checkbox"/> Safety pelvic belt		03972	£ 55
<input type="checkbox"/> Harness		04103	£ 109
<input type="checkbox"/> Lumbar support	(to insert in the backrest upholstery)	00013	£ 33
<input type="checkbox"/> Swing away lateral support	<input type="checkbox"/> Right	02097	£ 218
Not available with Posture backrest PHYSIO	<input type="checkbox"/> Left	02098	£ 218
<input type="checkbox"/> Spoke guards Black	(22"-24")	01043	£ 131
<input type="checkbox"/> Spoke guards Transparent	(22"-24")	00835	£ 131
<input type="checkbox"/> Spoke guards Progeo design	(22"-24")	00019	£ 153
<input type="checkbox"/> Spoke guards Dragon	(22"-24")	02106	£ 153




## Additional Kit

<input type="checkbox"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	£ 436
<input type="checkbox"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	£ 764

## PHYSIO - Posture backrest

Std. Backrest upholstery not provided

<input type="checkbox"/> Aluminum shell - Lateral retaining STANDARD 7,5cm	Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41 Height <input type="checkbox"/> 52	05308 05309	£ 884 £ 938
Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low			
<input type="checkbox"/> Aluminum shell - Lateral retaining DEEP 15cm	Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41 H25 not available for seat width 39-42-45-48 Height <input type="checkbox"/> 52	05310 05311	£ 938 £ 993
Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low			
<input type="checkbox"/> Carbon shell - Lateral retaining STANDARD 7,5cm	Height <input type="checkbox"/> 25 <input type="checkbox"/> 33 <input type="checkbox"/> 41	05312	£ 1015
Lateral retaining wing <input type="checkbox"/> A=high <input type="checkbox"/> B=low			




## COVER

☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




- ☐ SWING AWAY - Lateral support  
Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)




<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643
<input type="checkbox"/> Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)			

Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)



<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643

Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)



<input type="checkbox"/> Left	05646	<input type="checkbox"/> Right	05647	£ 175
<input type="checkbox"/> Left	05642	<input type="checkbox"/> Right	05643	£ 175
(with additional kit)				

Date .....	Order .....	English User name .....	Address .....
------------	-------------	----------------------------	---------------

☐ Headrest adjustable - Double ball mounting system



**Not available for Physio backrest's height 25 and with carbon shell**

- ☐ Large (34x14)  
☐ Small (28x14)



**05648 £ 393**

## Note

Italia Certified company normative UNI EN ISO 9001			Important notices page only if any optional have been chosen	
Pag. 5 / 5	Wheelchair : TEKNA ADVANCE - Rev. 2019.0	Prices in EURO - VAT excluded to be added to the price list	<input checked="" type="checkbox"/> standard	(supplied if no requests are indicated)
			<input type="checkbox"/> on request	(without extra charge)
			<input type="radio"/> optional	(cost to be added at the basic price)

## Enclosure to wheelchair order form Tekna Advance

### REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SMALL 	max 49	max 50	max 51	
MEDIUM 	max 53	max 54	max 55	max 57
The size of the fork will be defined only by Rehateam Progeo in accordance with wheelchair final configuration.				

2	Rear height			
	Rear wheels			
	22"	24"	25"	26"
H min	min. 36	min. 38	min. 39	min. 40
H max	max 44	max 48	max 49	max 51

3	Footplate option	Footplate distance
	Automatic Closure footplate	min. 36
	Automatic Closure footplate , with multifunction rotative system	min. 31
	One piece footplate - Flip-up!	min. 36
	One piece footplate - Flip-up!, with multifunction rotative system	min. 31
	One piece aluminium footplate - Flip-up!	min. 36
	Plastic footrest	min. 36
	One piece Carbon fibre footplate	min. 36
Calculated with the best configuration and with the footplates positioned on inside.		
<b>IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)</b>		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)