

TRILIFT

CLASSIC *Scooter & Wheelchair lift*



TRILIFT has hundreds of pre-fabricated customised docking brackets to meet your needs in the most affordable way.

Our patented system provides excellent scalability.

The TRILIFT integrated locking mechanism provides ease of use and ultimate security.

TRILIFT UK

NO BULKY PLATFORM



Reliable
Stabilised
High Performance
Automatically Locks & Secures

Unique
TRIangle Design

Fully Customised for Your Needs
The most ground clearance of any other outside lift




TRILIFT CLASSIC

HIGH PERFORMANCE



Tilt Adaptor gives you easy access to the rear door



Key switches provide extra security for lift operations.



The Red Triangle Device makes it easier to see and operate



Unique Triangular Design ensures better fit and ultimate stability

Benefits:

- No Platform
- No Straps
- Affordable
- Stabilised
- Easy to park
- More Ground Clearance
- Lower Profile Design
- Locked and Secured
- Fully Automatic
- Bolted Assembly
- Easy to Install
- Lightweight
- No Noise
- No Vibration
- Easy Operation
- 3 years warranty

NO MORE PLATFORM !!!

Contact your local dealer:
Mobility Supplies NZ
0800 000 605

- Transferable from vehicle to vehicle
- Heavy Duty Industrial Electric Motor
- Easily attaches to any vehicle with a tow bar including rental vehicles

Differences

TRILIFT CLASSIC

TRILIFT ULTRA-LITE

Lift Height	45" - 60"	42" - 55"
Lift Weight	60-80 lbs.	42-60 lbs.
Loading Capacity	up to 300 lbs.	up to 175 lbs.
Back Space	11"	10"
Cycle Load Time	55 sec	45 sec
Docking Bracket Size	9.5 "	7 "
Docking Bracket Weight	12-25 lbs.	8-15 lbs.

Options:

- 4 position Tilt Adaptor w/ spring load
- Coil-Cord Remote
- Hitch Pin Lock
- Lighted License Place Holder
- Soft Cover

TRILIFT is a Unique Lift & Lock Device
Works Easily with the PUSH OF A BUTTON



The NEXT GENERATION
in Mobility Transportation

TRILIFT -HD

*Scooter &
Wheelchair lift*



TRILIFT has hundreds of pre-fabricated customised docking brackets to meet your needs in the most affordable way.

Our patented system provides excellent scalability.

The TRILIFT integrated locking mechanism provides ease of use and ultimate security.

TRILIFT ^{UK}

450lb Loading Capacity



NO BULKY PLATFORM

Reliable

Stabilised

High Performance

Automatically Locks & Secures




**The most ground clearance of any other
outside lift**

TRILIFT HD

HIGH PERFORMANCE




Tilt Adaptor gives you easy access to the rear door



Key switches provide extra security for lift operations.



The Red Triangle Device makes it easier to see and operate



Unique Triangular Design ensures better fit and ultimate stability

Benefits:

- No Platform
- No Straps
- Affordable
- Stabilised
- Easy to park
- More Ground Clearance
- Lower Profile Design
- Locked and Secured
- Fully Automatic
- Bolted Assembly
- Easy to Install
- Lightweight
- No Noise
- No Vibration
- Easy Operation
- 3 years warranty

NO MORE PLATFORM !!!

Contact your local dealer:
Mobility Supplies NZ
0800 000 605

- Transferable from vehicle to vehicle
- Heavy Duty Industrial Electric Motor
- Easily attaches to any vehicle with a tow bar including rental vehicles

Differences	TRILIFT HD	TRILIFT CLASSIC
Lift Height	50" - 65"	45" - 60"
Lift Weight	65-95 lbs.	60-80 lbs.
Loading Capacity	up to 450 lbs.	up to 300 lbs.
Back Space	13"	11"
Cycle Load Time	55 sec	55 sec
Docking Bracket Size	9.5 "	9.5 "
Docking Bracket Weight	20-35 lbs.	12-25 lbs.

Options:

- 4 position Tilt Adaptor w/ spring load
- Coil-Cord Remote
- Hitch Pin Lock
- Lighted License Place Holder
- Soft Cover

TRILIFT is a Unique Lift & Lock Device
Works Easily with the PUSH OF A BUTTON



The NEXT GENERATION
in Mobility Transportation

The Lightest Lift on the market today!

TRILIFT ULTRA-LITE

*Accommodates Scooters Weighing
up to*

120lb



TRILIFT has hundreds of pre-fabricated customised docking brackets to meet your needs in the most affordable way.

Our patented system provides excellent scalability.

The TRILIFT integrated locking mechanism provides ease of use and ultimate security.

TRILIFT ^{UK}

**NO PLATFORM !
REVOLUTIONARY DESIGN**



Be More Independent!

Simple and Fully Automatic Operation.

LOCKS YOUR SCOOTER IN PLACE.

NO STRAPS!!!

**More Ground
Clearance**




**SIMPLE
& COMPACT**

TRILIFT ULTRA-LITE

THE LIGHTEST LIFT ON THE MARKET



Tilt Adaptor gives you easy access to the rear door



Benefits:

- No Platform
- No Straps
- Affordable
- Stabilised
- Easy to park
- More Ground Clearance
- Lower Profile Design
- Locked and Secured
- Fully Automatic
- Bolted Assembly
- Easy to Install
- Lightweight
- No Noise
- No Vibration
- Easy Operation
- 3 years warranty



The Red Triangle Device makes it easier to see and operate



Unique Triangular Design ensures better fit and ultimate stability

NO MORE HEAVY BULKY PLATFORM

Contact your local dealer:
Mobility Supplies NZ
0800 000 605

- Transferable from vehicle to vehicle
- Heavy Duty Industrial Electric Motor
- Easily attaches to any vehicle with a tow bar including rental vehicles

Differences

TRILIFT ULTRA-LITE

TRILIFT CLASSIC

Lift Height	42" - 55"	45" - 60"
Lift Weight	42-60 lbs.	60-80 lbs.
Loading Capacity	up to 175 lbs.	up to 300 lbs.
Back Space	10"	11"
Cycle Load Time	45 sec	55 sec
Docking Bracket Size	7 "	9.5 "
Docking Bracket Weight	8-15 lbs.	12-25 lbs.

Options:

- 4 position Tilt Adaptor w/ spring load
- Coil-Cord Remote
- Hitch Pin Lock
- Lighted License Place Holder
- Soft Cover

TRILIFT is a Unique Lift & Lock Device
Works Easily with the PUSH OF A BUTTON



The NEXT GENERATION
in Mobility Transportation